




CoverSeal™-AQ

AQUEOUS MOUNTING MEDIUM

 #FX2175-AQ 4OZ

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____

Formulated to be used after aqueous based stains, CoverSeal™-AQ is added as a pre-mounting medium (for tissue sections, cell spreads, cytopins or tissue imprints) following aqueous base stain prior to standard toluene or xylene mounting media. It has the advantage of tolerating water transferred from the stain procedure.



NON HAZARDOUS CHEMICAL

FIRST AID -

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical advice/attention.

IF ON SKIN: Remove/Take off immediately all contaminated clothing and wash it before reuse. Rinse skin with water/shower. If skin irritation occurs: get medical advice/attention. Call a POISON CENTER or doctor/physician if you feel unwell.

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Do NOT induce vomiting.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL -

Store locked up in a well ventilated place. Keep container tightly closed. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



CANCER
DIAGNOSTICS INC.

4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com **Made in USA**