



eosin TYPE w/PHLOXINE B



Eosin Type R w/Phloxine B is a classic formulation of the alcoholic eosin-Y counterstain with the addition of the Phloxine dye. It produces brilliant counterstaining with a more vivid red hue. Eosin Type R w/Phloxine B is ready to use and has been balanced to the proper pH for consistency and stability.

#EPHB 500ML

#EPHB-4 500ML 4/CS

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____



DANGER HAZARD STATEMENTS

Causes skin irritation. Causes serious eye irritation. Causes damage to organs (respiratory system). May cause damage to organs (respiratory system, central nervous system, liver, blood) through prolonged or repeated exposure. Flammable liquid and vapor.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye protection, face protection. Do not breathe dust, vapors. Do not eat, drink or smoke when using this product. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing If eye irritation persists: Get medical attention. If exposed or concerned: Call a doctor. Call a doctor if you feel unwell.

IF ON SKIN(or hair): Wash with plenty of water. Specific treatment (see first aid section on this label). If skin irritation or rash occurs: Get medical attention. Take off all contaminated clothing and wash it before reuse. Take off immediately all contaminated clothing. Rinse skin with water/shower.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.
STORAGE & DISPOSAL: Store locked up in a well ventilated place. Keep cool. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



CANCER
DIAGNOSTICS INC

4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com **Made in USA**