



eosin

TYPE V

routine counterstain



Eosin Type V is a classic formulation of the alcoholic eosin-Y counterstain. Eosin Type V produces brilliant counterstaining and differentiation of cell cytoplasm and other non-nuclear tissue components. Eosin Type V is ready to use and has been balanced to the proper pH for consistency and stability.

REF #ETV-4

500ML, Case of 4

LOT



For In Vitro Diagnostic Use

Date Received: _____ Date Opened: _____



DANGER HAZARD STATEMENTS

Causes skin irritation. Causes serious eye irritation. Causes damage to organs (respiratory system). May cause damage to organs (respiratory system, central nervous system, liver, blood) through prolonged or repeated exposure. Highly flammable liquid and vapor.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Do not breathe dust, vapors. Do not eat, drink, or smoke when using this product. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. If eye irritation persists: get medical attention. If exposed or concerned: call a doctor. Call a doctor if you feel unwell.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store locked up at room temperature in a well ventilated place. Keep container tightly closed. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



116 Page Point Circle Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**

POC284 Eosin Type V AWH-2741620 Rev.B 2023-12-26