



hematoxylin **T Y P E** progressive nuclear stain



The hybrid formulation of Hematoxylin Type R produces crisp nuclear detail and does not stain mucin. Hematoxylin Type R can be used in place of Gill's or Harris Hematoxylin formulations. This progressive stain allows users to adjust to their desired staining intensity by simply adjusting staining times. Hematoxylin Type R requires no daily filtering.

#HTR 500ML

#HTR-4 500ML 4/CS

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____



WARNING

HAZARD STATEMENTS

Causes skin irritation. Causes serious eye irritation. May cause damage to organs (kidneys) through prolonged or repeated exposure.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye protection, face protection. Do not breathe dust, vapors.

PRECAUTIONARY STATEMENTS - RESPONSE

IF ON SKIN: Wash with plenty of water. Specific treatment (see first aid section on this label). If skin irritation or rash occurs: Get medical attention. Take off all contaminated clothing and wash it before reuse.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing If eye irritation persists: Get medical attention. Call a doctor if you feel unwell.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL: Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com **Made in USA**