



hematoxylin **T Y P E** progressive nuclear stain



The unique formulation of Hematoxylin Type V produces crisp nuclear detail without background staining. Hematoxylin Type V can be used in place of Gill's Hematoxylin formulations. The progressive nature of this stain allows users to adjust to their desired staining intensity by simply adjusting staining times. Hematoxylin Type V requires no daily filtering

REF #HTV-4

500ML, Case of 4

LOT



For In Vitro Diagnostic Use

Date Received: _____ Date Opened: _____



WARNING

HAZARD STATEMENTS

Harmful if swallowed. May cause damage to kidneys through prolonged or repeated exposure (oral routes). Causes eye irritation.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash skin thoroughly after handling. Do not eat, drink or smoke when using this product. Do not breath mist/vapors/spray.

PRECAUTIONARY STATEMENTS - RESPONSE

IF SWALLOWED OR FEEL UNWELL: Call poison center or get medical advice/attention. Rinse mouth.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

STORAGE & DISPOSAL

Store locked up at room temperature. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



116 Page Point Circle Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**

POC283 Hematoxylin Type V AWH-2741623 Rev B 2023-12-26