



# ThinLayer

## Cell Preserve Solution Non-Gyn

REF #1NGYN-FF 32oz CS/4

LOT



IVD



59°F (15°C)  86°F (30°C)

### For In Vitro Diagnostics Use

Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_

Preservative solution for cell support during transport and preparation of slides. For use with liquid-based cytology method



#### **DANGER** **HAZARD STATEMENTS**

Causes skin irritation. Causes serious eye irritation. Causes damage to organs (respiratory system). May cause damage to organs (liver) through prolonged or repeated exposure. Flammable liquid and vapor.

#### **PRECAUTIONARY STATEMENTS - PREVENTION**

Wash body thoroughly after handling. Wear protective gloves. Wear eye protection, face protection. Do not breathe dust, vapors. Do not eat, drink or smoke when using this product. Keep away from heat sources and open flame. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical, ventilating, lighting and equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

#### **PRECAUTIONARY STATEMENTS - RESPONSE**

**IF ON SKIN (or hair):** Rinse/wash skin with plenty of water/shower. Specific treatment (see first aid section on this label). **If skin irritation or rash occurs:** Get medical attention. Take off immediately all contaminated clothing and wash it before reuse.

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. **If eye irritation persists, or exposed or concerned, or feel unwell:** Call a doctor/get medical attention.

**IN CASE OF FIRE:** Use water, dry chemical, CO<sub>2</sub> or foam to extinguish.

**STORAGE & DISPOSAL:** Store locked up, in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local regulations.

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



**CANCER**  
DIAGNOSTICS INC.

4300 Emperor Blvd - 400 Durham, NC 27703  
Tel: 1-877-846-5393 Fax: 1-877-817-1716  
[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)

MX1266 ThinLayer Cell Preserve Solution 32oz AWW-14105 Rev.B 2020-04-09