

MICROFLEX®

▲ BLACK DRAGON® ZERO



BLACK
NITRILE
NON-STERILE
AMBIDEXTROUS
TEXTURED FINGERTIPS
BEADED CUFF
FOR SINGLE USE ONLY
NON-LATEX

POWDER-FREE EXAMINATION GLOVES

If you are concerned about tactile sensitivity, consider Black Dragon® Zero gloves featuring textured fingertips. These gloves not only mask stains and provide better framing, they are also easy to put on and wear comfortably.

- Textured fingertips allow a secure grip on tools
- Black color masks stains and provides a backdrop for powders, solvents and other work elements



PRODUCT CODE		PACKAGING	
XS	BD-1000-NPF-XS	100	Gloves/Box
S	BD-1001-NPF-S	10	Boxes/Case
M	BD-1002-NPF-M		
L	BD-1003-NPF-L		
XL	BD-1004-NPF-XL		
XXL	BD-1005-NPF-XXL		

GLOVE PROFILE		ASTM D 6319	BLACK DRAGON® ZERO
Length	<i>mm (inches)</i>	min. 220 (8.7)	245 (9.6)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.07 (2.8)
Palm Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.09 (3.5)
Finger Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.14 (5.5)
Tensile Strength	<i>Before Aging</i>	min. 14	42
	<i>After Aging</i>	min. 14	14
Elasticity/Elongation	<i>Before Aging</i>	min. 500	500
	<i>After Aging</i>	min. 400	400
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED
INDUSTRIES