



10% Neutral Buffered Formalin

With Eosin Tint

REF #FX1003-T | Gallon (3.8L) CS/4

LOT



For In Vitro Diagnostic Use



Date Received: _____ Date Opened: _____

A general purpose fixative for all tissues to be examined by light microscopy. The addition of eosin to the 10% NBF formulation allows the user to more easily visualize small fragments of tissue during the transfer of the tissue from the fixative to the tissue cassette.



DANGER

HAZARD STATEMENTS

Harmful if inhaled. Causes skin irritation and serious eye damage. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction. May cause genetic defects. May cause cancer. May cause damage to organs (lungs, nose) through single, prolonged, or repeated exposure.

PRECAUTIONARY STATEMENTS - PREVENTION

Use only outdoors or in a well-ventilated area. Wash body thoroughly after handling. Wear protective gloves and clothing. Wear eye and face protection. Wear NIOSH approved respiratory protection. Contaminated work clothing must not be allowed out of the workplace. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Do not breathe dust, vapors. Do not eat, drink, or smoke when using this product.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. Immediately call a doctor.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off immediately all contaminated clothing and wash it before reuse.

IF INHALED: If breathing is difficult, remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: call a doctor. If exposed or concerned: get medical attention/call a doctor. Call a doctor if you feel unwell.

STORAGE & DISPOSAL

Store locked up at room temperature. Keep cool. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**

MX1977 NBF Eosin Tint Gallons, 4 AW-3524030 Rev.B 2022-03-14