



VaporPad™ Original

- Intended Use -

VaporPads™ are Formalin Neutralizer cutting pads that encapsulate and neutralize harmful vapors to greatly reduce exposure to Formaldehyde at the grossing station.

- General Information -

The VaporPad™ is constructed from a super absorbent nonwoven polyester composite. The surface is highly cut resistant due to the polyethylene dot laminate. Light blue contained work surface offers contrast for all biopsies. Cut resistant, durable and disposable.

- Specifications -

Width	up to 20	<i>inches</i>
Material	Polyester Fiber	
Weight	10.75	<i>oz/yd²</i>
Construction	Non-woven Composite	
Thickness	0.050 – 0.070 <i>inches</i>	<i>Inches</i>
Absorbency	10.0 mean	<i>gm/gm of fabric</i>
Laminate	Polyethylene Dot	<i>Top</i>
Color	Light Blue / White	<i>Top / Bottom</i>

- Packaging -

Catalog#	Size	Pack Qty
VP0810	8" x 10"	25
VP01114	11" x 14"	25
VP01620	16" x 20"	25
VP06014	16" x 60' Roll	Each
BLK0810	8" x 10" Bulk	400
BLK1114	11" x 14" Bulk	400

