



# VaporPad™ Ultra

## Cutting Pad

|                          |         |           |       |
|--------------------------|---------|-----------|-------|
| <input type="checkbox"/> | #UP0811 | 8" x 11"  | PK/25 |
| <input type="checkbox"/> | #UP1115 | 11" x 15" | PK/25 |

For In Vitro Diagnostics Use

Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_

VaporPad™ Ultra, Formalin Neutralizer Pads provide an effective way for lab personnel to guard against formaldehyde vapors and exposure by allowing formalin penetration into a highly absorbent core.

8" x 11" Pad will Absorb & Neutralize up to 120mL of Formalin  
11" x 15" Pad will Absorb & Neutralize up to 220mL of Formalin



### WARNING HAZARD STATEMENTS

Causes Skin Irritation. Causes serious eye irritation

### PRECAUTIONARY STATEMENTS - PREVENTION

Wash hands thoroughly after handling. Wear protective gloves / eye / face protection.

### PRECAUTIONARY STATEMENTS - RESPONSE

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical advice/attention.

**IF ON SKIN:** Rinse skin with plenty of water. Take off contaminated clothing and wash it before reuse.

### STORAGE & DISPOSAL

Keep bag tightly closed. Dispose of contents/container to an approved waste disposal plant.

**NOTICE:** Read Safety Data Sheet for additional information.

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



CANCER  
DIAGNOSTIC INC.

4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**