



Date Received: _____ Date Opened: _____

VaporPad™ Ultra, Formalin Neutralizer Pads provide an effective way for lab personnel to guard against formaldehyde vapors and exposure by allowing formalin penetration into a highly absorbent core.

8" x 11" Pad will Absorb & Neutralize up to 120mL of Formalin
11" x 15" Pad will Absorb & Neutralize up to 220mL of Formalin

VaporPad™ Ultra

Cutting Pad

<input type="checkbox"/>	#UP0811	8" x 11"	PK/25
<input type="checkbox"/>	#UP1115	11" x 15"	PK/25



WARNING HAZARD STATEMENTS

Causes Skin Irritation. Causes serious eye irritation

PRECAUTIONARY STATEMENTS - PREVENTION

Wash hands thoroughly after handling. Wear protective gloves / eye / face protection.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: get medical advice/attention.

IF ON SKIN: Rinse skin with plenty of water. Take off contaminated clothing and wash it before reuse.

STORAGE & DISPOSAL

Keep bag tightly closed. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300

For In Vitro Diagnostics Use



CANCER
DIAGNOSTIC INC.

4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com **Made in USA**