

Safety Data Sheet

Alpha Amylase



Section 1 Product Description

Product Name: Alpha Amylase
Recommended Use: Science, Laboratory or education applications Diastase
Synonyms:
Distributor: Cancer Diagnostics, Inc.
116 Page Point Circle Durham, NC 27703
1-877-846-5393
Chemtrec: 800-424-9300 (Transportation Spill Response 24 hours)

Section 2 Hazard Identification

Classification of the chemical in accordance with paragraph (d) of §1910.1200;

Not a dangerous substance according to GHS classification criteria. No known OSHA hazards.

GHS Classification:

Other Safety Precautions: May cause eye irritation.
May cause gastrointestinal discomfort.
May cause irritation to respiratory tract.
May cause irritation to skin.

Section 3 Composition / Information on Ingredients

<u>Chemical Name</u>	<u>CAS #</u>	<u>%</u>
Fungal alpha amylase	9013-01-8	100
Refrigerate upon arrival. {EMSFORM_03SDS_COMP_NOTE}		

Section 4 First Aid Measures

Emergency and First Aid Procedures

Inhalation: In case of accident by inhalation: remove casualty to fresh air and keep at rest.
Eyes: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
Skin Contact: After contact with skin, wash immediately with plenty of water.
Ingestion: If swallowed, do not induce vomiting: seek medical advice immediately and show this container or label.

Section 5 Firefighting Procedures

Extinguishing Media: Use media suitable to extinguish surrounding fire.
Fire Fighting Methods and Protection: Firefighters should wear full protective equipment and NIOSH approved self-contained breathing apparatus.
Fire and/or Explosion Hazards: N/A
Hazardous Combustion Products: Carbon dioxide, Carbon monoxide

Section 6 Spill or Leak Procedures

Steps to Take in Case Material Is Released or Spilled: No adverse health affects expected from the clean-up of spilled material. Ventilate the contaminated area.
Prevent the spread of any spill to minimize harm to human health and the environment if safe to do so. Wear complete and proper personal protective equipment following the recommendation of Section 8 at a minimum. Dike with suitable absorbent material like granulated clay. Gather and store in a sealed container pending a waste disposal evaluation. Avoid dusting. Contain the discharged material.

Safety Data Sheet

Section 7 Handling and Storage

Handling: Avoid creating and inhaling dust.
Storage: Keep container tightly closed in a cool, well-ventilated place.
Storage Code: Green - general chemical storage

Section 8 Protection Information

<u>Chemical Name</u>	<u>(TWA)</u>	<u>ACGIH</u>	<u>(STEL)</u>	<u>(TWA)</u>	<u>OSHA PEL</u>	<u>(STEL)</u>
No data available	N/A		N/A	N/A		N/A

Control Parameters

Engineering Measures: No exposure limits exist for the constituents of this product. General room ventilation might be required to maintain operator comfort under normal conditions of use.

Personal Protective Equipment (PPE): Lab coat, apron, eye wash, safety shower.

Respiratory Protection: No respiratory protection required under normal conditions of use.

Eye Protection: Wear chemical splash goggles when handling this product. Have an eye wash station available.

Skin Protection: Avoid skin contact by wearing chemically resistant gloves, an apron and other protective equipment depending upon conditions of use. Inspect gloves for chemical break-through and replace at regular intervals. Clean protective equipment regularly. Wash hands and other exposed areas with mild soap and water before eating, drinking, and when leaving work.

Gloves: Nitrile

Section 9 Physical Data

Formula: N/A	Vapor Pressure: N/A
Molecular Weight: N/A	Evaporation Rate (BuAc=1): N/A
Appearance: Off-white to tan Solid	Vapor Density (Air=1): >1
Odor: No data available	Specific Gravity: >1
Odor Threshold: No data available	Solubility in Water: Soluble
pH: No data available	Log Pow (calculated): No data available
Melting Point: No data available	Autoignition Temperature: No data available
Boiling Point: No data available	Decomposition Temperature: No data available
Flash Point: No data available	Viscosity: No data available
Flammable Limits in Air: N/A	Percent Volatile by Volume: N/A

Section 10 Reactivity Data

Reactivity: No data available
Chemical Stability: Stable under normal conditions.
Conditions to Avoid: None known.
Hazardous Polymerization: Will not occur

Section 11 Toxicity Data

Routes of Entry: Inhalation and skin contact.
Symptoms (Acute): Allergies
Delayed Effects: No data available

Acute Toxicity:

<u>Chemical Name</u>	<u>CAS Number</u>	<u>Oral LD50</u>	<u>Dermal LD50</u>	<u>Inhalation LC50</u>
No data available	9013-01-8	Not determined	Not determined	Not determined

Carcinogenicity:

<u>Chemical Name</u>	<u>CAS Number</u>	<u>IARC</u>	<u>NTP</u>	<u>OSHA</u>
No data available	9013-01-8	Not listed	Not listed	Not listed

Safety Data Sheet

Chronic Effects:

Mutagenicity: No evidence of a mutagenic effect.
Teratogenicity: No evidence of a teratogenic effect (birth defect).
Sensitization: No evidence of a sensitization effect.
Reproductive: No evidence of negative reproductive effects.
Target Organ Effects:
Acute: See Section 2
Chronic: N/A

Section 12

Ecological Data

Overview: This material is not expected to be harmful to the ecology.
Mobility: No data
Persistence: No data
Bioaccumulation: No data
Degradability: No data
Other Adverse Effects: No data

Chemical Name	CAS Number	Eco Toxicity
N/A	9013-01-8	

Section 13

Disposal Information

Disposal Methods: Dispose in accordance with all applicable Federal, State and Local regulations. Always contact a permitted waste disposer (TSD) to assure compliance.
Waste Disposal Code(s): Not Determined

Section 14

Transport Information

Ground - DOT Proper Shipping Name:
Not regulated for transport by US DOT.

Air - IATA Proper Shipping Name:
Not regulated for air transport by IATA.

Section 15

Regulatory Information

TSCA Status: All components in this product are on the TSCA Inventory.

Chemical Name	CAS Number	§ 313 Name	§ 304 RQ	CERCLA RQ	§ 302 TPQ	CAA 112(2) TQ
No data available	9013-01-8	No	No	No	No	No

California Prop 65:

No California Proposition 65 ingredients

Section 16

Additional Information

Safety Data Sheet

Revised: 2022-05-17

Replaces: 2019-01-11

The information provided in this (Material) Safety Data Sheet represents a compilation of data drawn directly from various sources available to us. Cancer Diagnostics, Inc. makes no representation or guarantee as to the suitability of this information to a particular application of the substance covered in the (Material) Safety Data Sheet.

Glossary

ACGIH	American Conference of Governmental Industrial Hygienists	NTP	National Toxicology Program
CAS	Chemical Abstract Service Number	OSHA	Occupational Safety and Health Administration
CERCLA	Comprehensive Environmental Response, Compensation, and Liability Act	PEL	Permissible Exposure Limit
DOT	U.S. Department of Transportation	ppm	Parts per million
IARC	International Agency for Research on Cancer	RCRA	Resource Conservation and Recovery Act
N/A	Not Available	SARA	Superfund Amendments and Reauthorization Act
		TLV	Threshold Limit Value
		TSCA	Toxic Substances Control Act
		IDLH	Immediately dangerous to life and health