



60below™ Freeze Spray



- **Non Flammable**
- Non-Ozone Depleting Rapidly
- Evaporating Zero Residue

Ozone Safe. For rapid freezing of fresh or paraffin embedded tissue.

60 Below™, Freeze Spray

Chills to -60F. Contains no CFC's or HCFC's. Straw attached to can to localize spray onto blocks and fresh tissue.

Item	Description	Qty
06212-NF	60 Below™	CS/12

10 oz (284g)

For in vitro diagnostic use only



Cancer Diagnostics Inc.

Contact Information

Visit us online

4300 Emperor Blvd-400
Durham, North Carolina 27703

Toll Free: (877) 846 5393
Fax: (877) 817 1716

cancerdiagnostics.com