

**Cutting forme**



Apply dyes to fresh or formalin fixed tissue. Only small amounts of dye are necessary for effective marking. Do not apply to living patients. Refer to SDS for more information.

**24 Hour Emergency Contact**  
**CHEMTREC USA and Canada 800-424-9300**

59°F (15°C) 86°F (30°C)

LOT

Cancer Diagnostics, Inc. 4300 Emperor Blvd - 400  
 Durham, NC 27703 USA

EC REP mdi Europa GmbH Langenhagener Str. 71  
 D-30855 Langenhagen Germany +49-511-3908 9530


**Warning** Causes eye irritation. Causes skin irritation. Wash skin thoroughly after handling. Wear eye protection/ face protection. Take off contaminated clothing and wash it before reuse. Store locked up. Dispose of contents/container in accordance with local regulations.

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**IF ON SKIN:** Wash with plenty of soap and water. If eye or skin irritation

IVD CE MADE IN USA

MX1171 Squeezer Box PK/15 CDi's® TISSUE MARKING DYE  
 117-AW-SQUEEZER\_BOX Rev.B 2019-09-09



**Features**

- ▶ **Optimal dye consistency**
- ▶ **Bright Dyes hold permanently through processing and staining**
- ▶ **Easy Application**
- ▶ **Packaging limits exposure to air**

- MD1001 BLACK
- MD1002 BLUE
- MD1003 GREEN
- MD1004 RED
- MD1005 YELLOW
- MD1006 ORANGE
- MD1007 VIOLET
- MD2000 (6 COLOR KIT; (3 OF EACH 6 COLORS (18), W/CAPS)
- MD3000 (7 COLOR KIT; (3 OF EACH 7 COLORS (21), W/CAPS)

cancerdiagnostics.com

