

Shake well before using

Apply dye to fresh or frozen fixed tissue. Only small amounts of dye are necessary for effective marking. Do not apply to living patients. Refer to SDS for more information.

RED

Warning: Causes skin irritation. Causes eye irritation. Wash skin thoroughly after handling. Wear eye and face protection. Take off contaminated clothing and wash it before reuse. Store Locked Up. Dispose of contents in accordance with local regulations.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **IF ON SKIN:** Wash with plenty of soap and water. If red or skin irritation persists, Get medical advice or attention.

 Cancer Diagnostics, Inc.
4200 Superior Blvd, #400, Durham NC 27705 USA
 and Europa GmbH Langenloper Str. 71 D-30855
Langenhagen Germany +49 511 3908 9330



1-800 (30°C)
800 (15°C)



24 Hour Emergency Contact: USA and Canada 800-424-9300
MX) 232 Applicator Series, Red AWW-17213 Rev.B 2020-02-05



CANCER
DIAGNOSTICS, INC.

CDI's®

Tissue Marking Dye
APPLICATOR SERIES™

REF 0726-4

.10 FL. OZ., (3ML), FINE LINE BRUSH CAP, PK/5

RED

cancerdiagnostics.com

MADE IN USA