

Sodium Thiosulfate Solution 5%


1. Identification

1.1 Product Identifier	Product Name: Sodium Thiosulfate Solution (5%) Product Number: SSC1131	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com

2. Hazards Identification

2.1 Classification: Regulation (EC) No 1272/2008 (CLP)	Not a hazardous mixture as classified by Regulation (EC) No 1272/2008 (CLP)
Classification: OSHA 29 CFR 1910.1200 HCS	Not a hazardous mixture as classified by OSHA 29 CFR 1910.1200 HCS (GHS)
Classification: Directive 1999/45/EC	Not a hazardous mixture as classified by Directive 1999/45/EC For full text of R- & S- phrases: see Section 15.

Sodium Thiosulfate Solution 5%

2.2 Label Elements	Pictogram: Not a hazardous mixture Signal word: Not a hazardous mixture. Hazard Statement(s): Not a hazardous mixture Precautionary Statements(s): Not a hazardous mixture		
NFPA Scale: 0-4			
HMIS (U.S.A.) Scale: 0-4	HEALTH	1	
	FLAMMABILITY	0	
	PHYSICAL HAZARD	0	
	PERSONAL PROTECTION	B	
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.		

3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Sodium Thiosulfate	7772-98-7	231-867-5	N/A	≤5

4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water. Get medical attention if irritation persists.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Sodium Thiosulfate Solution 5%

5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

6. Accidental Release Measures

<p>6.1 Personal precautions, protective equipment and emergency procedures Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.</p> <p>6.2 Environmental precautions Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.</p> <p>6.3 Methods and materials for containment and cleaning up Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.</p>
--

7. Handling and Storage

<p>7.1 Precautions for safe handling. Avoid contact with skin and eyes. Wash thoroughly after handling. Avoid breathing vapor.</p> <p>7.2 Conditions for safe storage, including any incompatibilities. Store in well ventilated area. Keep container tightly closed. Store at 15-30°C.</p> <p>7.3 Specific end use(s). See section 1.2</p>
--

8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: N/A
8.2 Exposure controls	<p>Personal Protective Equipment (PPE):</p> <p>Eye/Face protection. Safety glasses or goggles are required.</p> <p>Skin protection. Protective clothing is required.</p> <p>Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove.</p> <p>Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Sodium Thiosulfate Solution 5%

9. Physical and Chemical Properties

Physical State	Liquid
Color	Colorless
Odor	Odorless
Odor Threshold	Unknown
pH	Neutral
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Sodium nitrate, oxidizing agents, acids, alkalis.
10.6 Hazardous Decomposition Materials	Sodium oxide, Hydrogen sulfide.

11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. May be irritating to skin and mucous membranes. Serious Eye Damage/Irritation. May be irritating to the eyes. Respiratory or skin sensitization. No relevant data available.
---	---

Sodium Thiosulfate Solution 5%

	Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods. ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods.
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods. IATA: Not dangerous goods. IMDG: Not dangerous goods. ADR/RID: Not dangerous goods.
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Sodium Thiosulfate Solution 5%

15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R36/38 EN: Irritating to eyes and skin. DA: Irriterer øjnene og huden DE: Reizt die Augen und die Haut. EL: Ερεθίζει τα μάτια και το δέρμα. ES: Irrita los ojos y la piel. FR: Irritant pour les yeux et la peau. IT: Irritante per gli occhi e la pelle. NL: Irriterend voor de ogen en de huid. PT: Irritante para os olhos e pele. SV: Irriterar ögonen och huden.
S phrases	S2 EN: Keep out of the reach of children. DA: Opbevares utilgængeligt for børn. DE: Darf nicht in die Hände von Kindern gelangen. EL: Μακριά από παιδιά. ES: Manténgase fuera del alcance de los niños. FR: Conserver hors de la portee des enfants. IT: Conservare fuori della portata dei bambini. NL: Buiten bereik van kinderen bewaren. PT: Manter fora do alcance das crianças. SV: Förvaras oåtkomligt för barn. S24/25 EN: Avoid contact with skin and eyes. DA: Undgå kontakt med huden og øjnene. DE: Berührung mit den Augen und der Haut vermeiden. EL: Αποφεύγετε επαφή με το δέρμα και με τα μάτια. ES: Evítese el contacto con los ojos y la piel. FR: Éviter le contact avec la peau et les yeux. IT: Evitare il contatto con gli occhi e con la pelle. NL: Aanraking met de ogen en de huid vermijden. PT: Evitar o contacto com a pele e os olhos. SV: Undvik kontakt med huden och ögonen.

16. Other Information

Revision Date: 2019-01-11

Notice to Reader:


To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Silver Nitrate Solution 10%




1. Identification

1.1 Product Identifier	Product Name: Silver Nitrate Solution (10%) Product Number: SSC1121	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com

2. Hazards Identification

2.1 GHS Classification	Skin irritation (Category 2) – H315 Eye irritation (Category 2A) – H319 Acute aquatic toxicity (Category 1) – H400 Chronic aquatic toxicity (Category 1) – H410	
Classification: Directive 1999/45/EC	 <p>Xi EN: Irritant DA: Lokalirriterende DE: Reizend EL: Ερεθιστικό ES: Irritante FR: Irritant IT: Irritante NL: Irriterend PT: Irritante SV: Irriterande</p>	

Silver Nitrate Solution 10%

	 <p>N: EN: Dangerous for the environment DA: Miljøfarlig DE: Umweltgefährlich EL: Επικίνδυνο για το περιβάλλον ES: Peligroso para el medio ambiente FR: Dangereux pour l'environnement IT: Pericoloso per l'ambiente NL: Milieugevaarlijk PT: Perigoso para o ambiente SV: Miljöfarlig</p> <p>For full text of R- & S- phrases: see Section 15.</p>								
2.2 Label Elements	<p>Pictogram(s):</p>  <p>Signal word: Warning</p> <p>Hazard statement(s): H315 – Causes skin irritation. H319 – Causes serious eye irritation H410 – Very toxic to aquatic life with long lasting effects.</p> <p>Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 + - IF IN EYES: Rinse cautiously with water for several minutes. P337 + P 313 - If eye irritation persists get medical advice/attention.</p>								
NFPA Scale: 0-4									
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr> <td>HEALTH</td> <td>1</td> </tr> <tr> <td>FLAMMABILITY</td> <td>0</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>B</td> </tr> </table>	HEALTH	1	FLAMMABILITY	0	PHYSICAL HAZARD	0	PERSONAL PROTECTION	B
HEALTH	1								
FLAMMABILITY	0								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	B								
2.3 Other Hazards	<p>PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.</p>								





Silver Nitrate Solution 10%

3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Silver Nitrate	7761-88-8	231-853-9	 Danger. 2. H272 May intensify fire; oxidizer.  Warning. 4 H302 Harmful if swallowed.  Danger. 1B H314 Causes severe skin burns and eye damage.  Warning. H410 Very toxic to aquatic life with long lasting effects.	≤ 10

4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water. Get immediate medical attention if irritation persists.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	None known
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Silver Nitrate Solution 10%

6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 2-8°C.

7.3 Specific end use(s).

See section 1.2

8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: N/A
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

9. Physical and Chemical Properties

Physical State	Liquid
Color	Transparent, colorless
Odor	Unknown
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown

Silver Nitrate Solution 10%

Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Oxidizers, alkalis, antimony salts, arsenates, carbonates and phosphates.
10.6 Hazardous Decomposition Materials	Carbon oxides, nitrogen oxides, oxides of silver.

11. Toxicological Information

11.1 Information on Toxicological Effects.	<p>Acute Toxicity. No relevant data available</p> <p>Skin Corrosion/Irritation. Irritating to skin and mucous membranes.</p> <p>Serious Eye Damage/Irritation. Seriously irritating to eye.</p> <p>Respiratory or skin sensitization. No relevant data available.</p> <p>Germ Cell Mutagenicity. No relevant data available.</p>
Carcinogenicity.	<p>International Agency for Research on Cancer (IARC). None of the components are listed.</p> <p>National Toxicology Program (NTP). None of the components are listed.</p>

Silver Nitrate Solution 10%

12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: 3082 ADR/RID: 3082 IMDG: 3082 IATA: 3082
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Environmentally hazardous substance, liquid, n.o.s.(Silver nitrate) IMDG: Environmentally hazardous substance, liquid, n.o.s.(Silver nitrate) IATA: Environmentally hazardous substance, liquid, n.o.s.(Silver nitrate) DOT: Environmentally hazardous substance, liquid, n.o.s.(Silver nitrate)
14.3 Transport Hazard Class(es)	DOT: 9 IATA: 9 IMDG: 9 ADR/RID: 9
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: III IMDG: III IATA: III DOT: III
14.5 Environmental Hazards	Marine Pollutant: Yes
14.6 Special Precautions for User	Not applicable.

Silver Nitrate Solution 10%

15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	<p>R36/38 EN: Irritating to eyes and skin. DA: Irriterer øjnene og huden DE: Reizt die Augen und die Haut. EL: Ερεθίζει τα μάτια και το δέρμα. ES: Irrita los ojos y la piel. FR: Irritant pour les yeux et la peau. IT: Irritante per gli occhi e la pelle. NL: Irriterend voor de ogen en de huid. PT: Irritante para os olhos e pele. SV: Irriterar ögonen och huden.</p> <p>R50/53 EN: Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment. DA: Meget giftig for organismer, der lever i vand; kan forårsage uønskede langtidsvirkninger i vandmiljøet. DE: Sehr giftig für Wasserorganismen, kann in Gewässern längerfristig schädliche Wirkungen haben. EL: Πολύ τοξικό για τους υδρόβιους οργανισμούς, μπορεί να προκαλέσει μακροχρόνιες δυσμενείς επιπτώσεις στο υδάτινο περιβάλλον. ES: Muy tóxico para los organismos acuáticos, puede provocar a largo plazo efectos negativos en el medio ambiente acuático. FR: Très toxique pour les organismes aquatiques, peut entraîner des effets néfastes a long terme pour l'environnement aquatique. IT: Altamente tossico per gli organismi acquatici, può provocare a lungo termine effetti negativi per l'ambiente acquatico. NL: Zeer vergiftig voor in het water levende organismen; kan in het aquatisch milieu op lange termijn schadelijke effecten veroorzaken. PT: Muito tóxico para os organismos aquáticos, podendo causar efeitos nefastos a longo prazo no ambiente aquático. SV: Mycket giftigt för vattenlevande organismer, kan orsaka skadliga långtidseffekter i vattenmiljön.</p>
S phrases	<p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspolien en deskundig medisch advies inwinnen.</p>

Silver Nitrate Solution 10%

	<p>PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista.</p> <p>SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S45</p> <p>EN: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).</p> <p>DA: Ved ulykkestilfælde eller ved ildebefindende er omgående lægebehandling nødvendig; vis etiketten, hvis det er muligt.</p> <p>DE: Bei Unfall oder Unwohlsein sofort Arzt zuziehen (wenn möglich, dieses Etikett vorzeigen).</p> <p>EL: Σε περίπτωση ατυχήματος ή αν αισθανθείτε αδιαθεσία ζητήστε αμέσως ιατρική συμβουλή (δείξτε την ετικέτα αν είναι δυνατό).</p> <p>ES: En caso de accidente o malestar, acúdase inmediatamente al médico (si es posible, muéstresela la etiqueta).</p> <p>FR: En cas d'accident ou de malaise, consulter immédiatement un medecin (si possible lui montrer l'etiquette).</p> <p>IT: In caso di incidente o di malessere consultare immediatamente il medico (se possibile, mostrargli l'etichetta).</p> <p>NL: Bij een ongeval of indien men zich onwel voelt, onmiddellijk een arts raadplegen (indien mogelijk hem dit etiket tonen).</p> <p>PT: Em caso de acidente ou de indisposição, consultar imediatamente o médico (se possível mostrar-lhe o rótulo).</p> <p>SV: Vid olycksfall, illamående eller annan påverkan, kontakta omedelbart läkare. Visa om möjligt etiketten.</p> <p>S60</p> <p>EN: This material and its container must be disposed of as hazardous waste.</p> <p>DA: Undgå udledning til miljøet. Se særlig vejledning/leverandørbrugsanvisning.</p> <p>DE: Dieses Produkt und sein Behälter sind als gefährlicher Abfall zu entsorgen.</p> <p>EL: Αποφύγετε την ελευθέρωσή του στο περιβάλλον. Αναφερθείτε σε ειδικές οδηγίες/.ελτίο δεδομένων ασφαλείας.</p> <p>ES: Elimínense el producto y su recipiente como residuos peligrosos.</p> <p>FR: Eliminer le produit et son recipient comme un dechet dangereux.</p> <p>IT: Questo materiale e il suo contenitore devono essere smaltiti come rifiuti pericolosi.</p> <p>NL: Deze stof en de verpakking als gevaarlijk afval afvoeren.</p> <p>PT: Este produto e o seu recipiente devem ser eliminados como resíduos perigosos.</p> <p>SV: Detta material och dess behållare skall tas om hand som farligt avfall.</p> <p>S61</p> <p>EN: Avoid release to the environment. Refer to special instructions/Safety data sheets.</p> <p>DA: Undgå udledning til miljøet. Se særlig vejledning/leverandørbrugsanvisning.</p> <p>DE: Freisetzung in die Umwelt vermeiden. Besondere Anweisungen einholen/Sicherheitsdatenblatt zu Rate ziehen.</p> <p>EL: Αποφύγετε την ελευθέρωσή του στο περιβάλλον. Αναφερθείτε σε ειδικές οδηγίες/.ελτίο δεδομένων ασφαλείας.</p> <p>ES: Evítese su liberación al medio ambiente. Recábense instrucciones específicas de la ficha de datos de seguridad.</p> <p>FR: Eviter le rejet dans l'environnement. Consulter les instructions speciales/la fiche de donnees de securite.</p> <p>IT: Non disperdere nell'ambiente. Riferirsi alle istruzioni speciali/ schede informative in materia di sicurezza.</p> <p>NL: Voorkom lozing in het milieu. Vraag om speciale instructies/veiligheidskaart.</p> <p>PT: Evitar a libertação para o ambiente. Obter instruções específicas/fichas de segurança.</p> <p>SV: Undvik utsläpp till miljön. Läs särskilda instruktioner/varuinformationsblad.</p>
--	--

16. Other Information

Revision Date: 2019-01-11

Notice to Reader:


To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Nuclear Fast Red Solution



1. Identification

1.1 Product Identifier	Product Name: Nuclear Fast Red Solution (Enhanced Stability) Product Number: SSC1099	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com

2. Hazards Identification

2.1 GHS Classification	Skin irritation (Category 2) – H315 Serious eye damage (Category 1) – H318 Chronic aquatic toxicity (Category 3) – H412	
Classification: Directive 1999/45/EC	 <p>Xi, R36, S26 EN: Irritant DA: Lokalirriterende DE: Reizend EL: Ερεθιστικό ES: Irritante FR: Irritant IT: Irritante NL: Irriterend PT: Irritante SV: Irriterande</p> <p>For full text of R- & S- phrases: see Section 15.</p>	

Nuclear Fast Red Solution




2.2 Label Elements	Pictogram(s):  Signal word: Danger Hazard statement(s): H303 – May be harmful if swallowed. H318 – Causes serious eye damage H412 – Harmful to aquatic life with long-lasting effects Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 + P315 - IF IN EYES: Rinse cautiously with water for several minutes. Get immediate medical advice/attention. P273 – Avoid release to the environment.			
NFPA Scale: 0-4				
HMIS (U.S.A.) Scale: 0-4	HEALTH	2		
	FLAMMABILITY	0		
	PHYSICAL HAZARD	0		
	PERSONAL PROTECTION	B		
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.			

3. Composition and Information on Ingredients




3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Water	7732-18-5	231-791-2	n/a	80-100
Aluminum Sulfate Hydrate	7784-31-8	233-135-0	 Warning. 2 H315 Causes skin irritation. 3 H335 May cause respiratory irritation.  Danger. 2 H318 Causes serious eye damage.  Warning. 2 H411 Toxic to aquatic life with long lasting effects.	≤ 5

Nuclear Fast Red Solution

Nuclear Fast Red	6409-77-4	229-088-0	 Warning. 2 H315 Causes skin irritation. 2 H319 Causes serious eye irritation. 3 H335 May cause respiratory irritation.	≤ 0.2
Sodium Azide	26628-22-8	247-852-1	 Danger 1, 2 H300 Fatal if swallowed., H310 Fatal in contact with skin.  Warning. 1. H410 Very toxic to aquatic life with long lasting effects.	≤ 0.1

4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	Sulphur oxides, Aluminum oxide
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Nuclear Fast Red Solution

7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.
 Wash thoroughly after handling.
 Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.
 Keep container tightly closed.
 Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Sodium Azide: ACGIH TLV: 0.29/mg/m ³ ceiling NIOSH REL-C: 0.3mg/m ³ skin
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

9. Physical and Chemical Properties

Physical State	Liquid
Color	Red
Odor	Odorless
Odor Threshold	Unknown
pH	3.2 ± 0.4
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown

Nuclear Fast Red Solution

Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Strong oxidizing agents, metals.
10.6 Hazardous Decomposition Materials	Sulphur oxides, Aluminum oxide.

11. Toxicological Information

11.1 Information on Toxicological Effects.	<p>Acute Toxicity. No relevant data available</p> <p>Skin Corrosion/Irritation. Irritating to skin and mucous membranes.</p> <p>Serious Eye Damage/Irritation. Causes serious damage to eye.</p> <p>Respiratory or skin sensitization. No relevant data available.</p> <p>Germ Cell Mutagenicity. No relevant data available.</p>
Carcinogenicity.	<p>International Agency for Research on Cancer (IARC). None of the components are listed.</p> <p>National Toxicology Program (NTP). None of the components are listed.</p>

12. Ecological Information

12.1 Toxicity	<p>Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.</p> <p>Other: Known to be harmful to aquatic life with long-lasting effects.</p>
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.

Nuclear Fast Red Solution

12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods. ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods.
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods. IATA: Not dangerous goods. IMDG: Not dangerous goods. ADR/RID: Not dangerous goods.
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.
14.7 Transport in Bulk According to Annex II of MARPOL73/78 and the IBC Code.	This product is not intended to be transported in bulk as defined in Annex II of MARPOL73/78 and the IBC Code.

15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	This mixture does not contain a chemical known to the state of California to cause cancer.
Right To Know Components	None of the components in this mixture are listed.

Nuclear Fast Red Solution

15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	<p>R38 EN: Irritating to skin. DA: Irriterer huden. DE: Reizt die Haut. EL: Ερεθίζει το δέρμα. ES: Irrita la piel. FR: Irritant pour la peau. IT: Irritante per la pelle. NL: Irriterend voor de huid. PT: Irriterend voor de huid. SV: Irriterar huden.</p> <p>R41 EN: Risk of serious damage to eyes. DA: Risiko for alvorlig øjenskade. DE: Gefahr ernster Augenschäden. EL: Κίνδυνος σοβαρών οφθαλμικών βλαβών. ES: Riesgo de lesiones oculares graves. FR: Risque de lésions oculaires graves. IT: Rischio di gravi lesioni oculari. NL: Gevaar voor ernstig oogletsel. PT: Risco de lesões oculares graves. SV: Risk för allvarliga ögonskador.</p> <p>R52 EN: Harmful to aquatic organisms. DA: Skadelig for organismer, der lever i vand. DE: Schädlich für Wasserorganismen. EL: Επιβλαβές για τους υδρόβιους οργανισμούς. ES: Nocivo para los organismos acuáticos. FR: Nocif pour les organismes aquatiques. IT: Nocivo per gli organismi acquatici. NL: Schadelijk voor in het water levende organismen. PT: Nocivo para os organismos aquáticos. SV: Skadligt för vattenlevande organismer.</p>
S phrases	<p>S36/37/39 EN: Wear suitable protective clothing, gloves and eye/face protection. DA: Brug særligt arbejdstøj, egnede beskyttelseshandsker og -briller/ansigtsskærm. DE: Bei der Arbeit geeignete Schutzkleidung, Schutzhandschuhe und Schutzbrille/Gesichtsschutz tragen. EL: Φοράτε κατάλληλη προστατευτική ενδυμασία, γάντια και συσκευή προστασίας ματιών/προσώπου ES: Úsense indumentaria y guantes adecuados y protección para los ojos/la cara. FR: Porter un vêtement de protection approprié, des gants et un appareil de protection des yeux/du visage. IT: Usare indumenti protettivi e guanti adatti e proteggersi gli occhi/la faccia. NL: Draag geschikte beschermende kleding, handschoenen en een beschermingsmiddel voor de ogen/het gezicht. PT: Usar vestuário de protecção, luvas e equipamento protector para os olhos /face adequados. SV: Använd lämpliga skyddskläder, skyddshandskar samt skyddsglasögon eller ansiktsskydd.</p> <p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico.</p>

Nuclear Fast Red Solution

	<p>FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste.</p> <p>IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico.</p> <p>NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen.</p> <p>PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista.</p> <p>SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p>
--	--

16. Other Information

Revision Date: 01/11/2019

Notice to Reader:



To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Gold Chloride Solution (0.2%) Product Number: SSC1080	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 Classification: Regulation (EC) No 1272/2008 (CLP)	Skin sensitization (Category 1) – H317
Classification: OSHA 29 CFR 1910.1200 HCS	Skin sensitization (Category 1) – H317
Classification: Directive 1999/45/EC	R42/43, S26, S36/37/39, S45 EN: May cause sensitization by inhalation and skin contact. DA: Kan give overfølsomhed ved indånding og ved kontakt med huden. DE: Sensibilisierung durch Einatmen und Hautkontakt möglich. EL: Μπορεί να προκαλέσει ευαισθητοποίηση όταν εισπνέεται και σε επαφή με το δέρμα. ES: Posibilidad de sensibilización por inhalación y por contacto con la piel. FR: Peut entraîner une sensibilisation par inhalation et par contact avec la peau. IT: Può provocare sensibilizzazione per inalazione e contatto con la pelle. NL: Kan overgevoeligheid veroorzaken bij inademing of contact met de huid. PT: Pode causar sensibilização por inalação e em contacto com a pele. SV: Kan ge allergi vid inandning och hudkontakt. For full text of R- & S- phrases: see Section 15.


2.2 Label Elements	Pictogram(s):  Signal word: Warning Hazard statement(s): H317 – May cause an allergic skin reaction Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes.		
NFPA Scale: 0-4			
HMIS (U.S.A.) Scale: 0-4	HEALTH	1	
	FLAMMABILITY	0	
	PHYSICAL HAZARD	0	
	PERSONAL PROTECTION	B	
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.		

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Gold Chloride	16961-25-4	240-948-4		≤ 0.2

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water. Get immediate medical attention if irritation continues.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	Decomposes on heating to chlorine gas, hydrogen chloride and metallic gold.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures
6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage
7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 2-8°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: N/A
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor.

	Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Transparent yellow
Odor	Odorless
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Ammonia or ammonia salts. Strong bases.
10.6 Hazardous Decomposition Materials	Hydrogen chloride.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Irritating to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods

14.3 Transport Hazard Class(es)	DOT: Not dangerous goods IATA: Not dangerous goods IMDG: Not dangerous goods ADR/RID: Not dangerous goods
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R42/43 EN: May cause sensitization by inhalation and skin contact. DA: Kan give overfølsomhed ved indånding og ved kontakt med huden. DE: Sensibilisierung durch Einatmen und Hautkontakt möglich. EL: Μπορεί να προκαλέσει ευαισθητοποίηση όταν εισπνέεται και σε επαφή με το δέρμα. ES: Posibilidad de sensibilización por inhalación y por contacto con la piel. FR: Peut entraîner une sensibilisation par inhalation et par contact avec la peau. IT: Può provocare sensibilizzazione per inalazione e contatto con la pelle. NL: Kan overgevoeligheid veroorzaken bij inademing of contact met de huid. PT: Pode causar sensibilização por inalação e em contacto com a pele. SV: Kan ge allergi vid inandning och hudkontakt.
S phrases	S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen. PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista. SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.

	<p>S36/37/39 EN: Wear suitable protective clothing, gloves and eye/face protection. DA: Brug særligt arbejdstøj, egnede beskyttelseshandsker og -briller/ansigtsskærm. DE: Bei der Arbeit geeignete Schutzkleidung, Schutzhandschuhe und Schutzbrille/Gesichtsschutz tragen. EL: Φοράτε κατάλληλη προστατευτική ενδυμασία, γάντια και συσκευή προστασίας ματιών/προσώπου ES: Úsense indumentaria y guantes adecuados y protección para los ojos/la cara. FR: Porter un vêtement de protection approprié, des gants et un appareil de protection des yeux/du visage. IT: Usare indumenti protettivi e guanti adatti e proteggersi gli occhi/la faccia. NL: Draag geschikte beschermende kleding, handschoenen en een beschermingsmiddel voor de ogen/het gezicht. PT: Usar vestuário de protecção, luvas e equipamento protector para os olhos /face adequados. SV: Använd lämpliga skyddskläder, skyddshandskar samt skyddsglasögon eller ansiktsskydd.</p> <p>S45 EN: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). DA: Ved ulykkestilfælde eller ved ildebefindende er omgående lægebehandling nødvendig; vis etiketten, hvis det er muligt. DE: Bei Unfall oder Unwohlsein sofort Arzt zuziehen (wenn möglich, dieses Etikett vorzeigen). EL: Σε περίπτωση ατυχήματος ή αν αισθανθείτε αδιαθεσία ζητήστε αμέσως ιατρική συμβουλή (δείξτε την επκέτα αν είναι δυνατό). ES: En caso de accidente o malestar, acúdase inmediatamente al médico (si es posible, muéstresele la etiqueta). FR: En cas d'accident ou de malaise, consulter immédiatement un medecin (si possible lui montrer l'etiquette). IT: In caso di incidente o di malessere consultare immediatamente il medico (se possibile, mostrargli l'etichetta). NL: Bij een ongeval of indien men zich onwel voelt, onmiddellijk een arts raadplegen (indien mogelijk hem dit etiket tonen). PT: Em caso de acidente ou de indisposição, consultar imediatamente o médico (se possível mostrar-lhe o rótulo). SV: Vid olycksfall, illamående eller annan påverkan, kontakta omedelbart läkare. Visa om möjligt etiketten.</p>
--	---

Section 16. Other Information

Revision Date: 2019-01-11

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.