

Date Received: _____ Date Opened: _____



Scan For Protocol

Movat Pentachrome Stain Kit

REF #SSI037-VO

This Movat Pentachrome Stain Kit (Modified Russel-Movat) is intended for use in the histological demonstration of collagen, elastin, muscle, and mucin in FFPE sections. This procedure has historically been particularly useful when studying the heart, blood vessels and various vascular diseases. A trichrome-type of staining using Biebrich Scarlet/Acid, Phosphotungstic Acid, and a proprietary Yellow Stain Solution, is then employed giving red muscle and yellow collagen.

For In Vitro Diagnostic Use



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

HC2000 Movat Pentachrome VO (L) (W) (M) (S) Rev. B 2021.10.21