

Date Received: _____ Date Opened: _____



Scan For Protocol

Orcein Stain Kit

(For Hepatitis B and Elastic Fibers)

REF #SSI039-VO

The Orcein Stain may be used in histology procedures for the visualization of Hepatitis B surface Antigen (HBsAg), elastic fibers, and copper associated proteins. HBsAg appears as irregular shaped aggregates in the cytoplasmic region of the cells. This reagent may be used on formalin-fixed, paraffin- embedded or frozen sections.

For In Vitro Diagnostic Use



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

PK2021-Orcein VO-01 (4/1/2021) 1.0 Rev B 2021.10.21