

Date Received: _____ Date Opened: _____



Scan For Protocol

Picro-Sirius Red Stain Kit (For Cardiac Muscle)

REF #SSI038-CM-VO

The Picro-Sirius Red Stain Kit (For Cardiac Muscle) is intended for use in the histological visualization of thin septa and collagen fibers. This modification of the Picro-Sirius Red stain eliminates the yellow cytoplasmic staining that can obscure thin collagenous septa. Using this procedure easily allows viewing of collagenous septa as thin as 0.2-0.5 microns. The PSR stain may be viewed using standard light microscopy or polarized light resulting in birefringence of the collagen fibers to distinguish between type I and type III.

For In Vitro Diagnostic Use



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

PSR238-Picro-Sirius-Red-Cardiac-Muscle-VO-1K-ASV-08/07/11 Rev. 8/2012-05/13