



GMS Fungus

Stain Kit

#SS1008-125 125ML

#SS1008-250 250ML

#SS1008-500 500ML

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____

The GMS (Grocott's Methenamine Silver) Fungus Stain Kit is a modified silver stain used to distinguish fungi, basement membranes and opportunistic organisms. Fungi and *Pneumocystis carinii* stain gray to black and mucin appears taupe to dark gray. These are present on a green background counterstained by Light Green.

Kit Includes:

BOX 1

- 1 x Chromic Acid, 5%, 250mL
- 1 x Sodium Borate, 5%, 250mL
- 1 x Sodium Thiosulfate, 2%, 250mL
- 1 x Light Green Stock Solution, 250mL
- 1 x Sodium Bisulfite, 1%, 250mL

BOX 2 - KEEP REFRIGERATED

- 1 x Gold Chloride, 0.1%, 250mL
- 1 x Methenamine-Silver Stock Solution, 250mL

STORE AT ROOM TEMPERATURE



4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716
www.cancerdiagnostics.com **Made in USA**



GMS Fungus

Stain Kit

#SS1008-125 125ML

#SS1008-250 250ML

#SS1008-500 500ML

For In Vitro Diagnostics Use

Date Received: _____ Date Opened: _____

The GMS (Grocott's Methenamine Silver) Fungus Stain Kit is a modified silver stain used to distinguish fungi, basement membranes and opportunistic organisms. Fungi and *Pneumocystis carinii* stain gray to black and mucin appears taupe to dark gray. These are present on a green background counterstained by Light Green.

Kit Includes:

BOX 1

- 1 x Chromic Acid, 5%, 250mL
- 1 x Sodium Borate, 5%, 250mL
- 1 x Sodium Thiosulfate, 2%, 250mL
- 1 x Light Green Stock Solution, 250mL
- 1 x Sodium Bisulfite, 1%, 250mL

BOX 2 - KEEP REFRIGERATED

- 1 x Gold Chloride, 0.1%, 250mL
- 1 x Methenamine-Silver Stock Solution, 250mL

KEEP REFRIGERATED (2°-8° C)



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**