





Section 1. Identification

1.1 Product Identifier	Product Name: Bouin's Fluid Product Number: SSC1053	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification










2.1 GHS Classification	Acute toxicity, Oral (Category 4), H302 Skin corrosion (Category 1A), H314 Serious eye damage (Category 1), H314 Carcinogenicity (Category 2), H351 Specific target organ toxicity - single exposure (Category 3), Respiratory system, H335
Classification: Directive 1999/45/EC	 Xn, R4, R20/21/22, R45; S26, S45, S51, S60 EN: Harmful DA: Sundhedsskadelig DE: Gesundheitsschädlich EL: Επιβλαβές ES: Nocivo FR: Nocif IT: Nocivo NL: Schadelijk

	<p>PT: Nocivo SV: Hälsoskadlig</p>  <p>C EN: Corrosive DA: Ætsende DE: Ätzend EL: Διαβρωτικό ES: Corrosivo FR: Corrosif IT: Corrosivo NL: Bijtend PT: Corrosivo SV: Frätande</p> <p>For full text of R- & S- phrases: see Section 15.</p>								
2.2 Label Elements	<p>Pictogram(s):</p>  <p>Signal word: Warning</p> <p>Hazard statement(s): H302 – Harmful if swallowed. H314 – Causes severe skin burns and eye damage H351 – Suspected of causing cancer H335 – May cause respiratory irritation</p> <p>Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes. P301 + P310 + P305 - IF SWALLOWED / IF IN EYES: Immediately call a POISON CENTER or doctor/ physician</p>								
NFPA Scale: 0-4									
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr> <td>HEALTH</td> <td>3</td> </tr> <tr> <td>FLAMMABILITY</td> <td>1</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>H</td> </tr> </table>	HEALTH	3	FLAMMABILITY	1	PHYSICAL HAZARD	0	PERSONAL PROTECTION	H
HEALTH	3								
FLAMMABILITY	1								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	H								
2.3 Other Hazards	<p>PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.</p>								

Section 3. Composition and Information on Ingredients
3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Formaldehyde	50-00-0	200-001-8	 Warning. 2 H315 Causes skin irritation.  Warning. 1. H400 Very toxic to aquatic life . 1. H410 Very toxic to aquatic life with long lasting effects.	≤ 8
Acetic Acid	64-19-7	200-580-7	 Warning. 3, H226 Flammable liquid and vapour.  Danger. 1A H314 Causes severe skin burns and eye damage.	≤ 5
Methyl Alcohol	67-56-1	200-659-6	 Danger. 3 H225 Highly flammable liquid and vapour.  Danger. 3. H301 Toxic if swallowed. 3 H311 Toxic in contact with skin. 3 H331 Toxic if inhaled.  Danger. 1 H370 Causes damage to organs.	≤ 2
Picric Acid	88-89-1	201-865-9	 Danger. Div 1.1 H201 Explosive; mass explosion hazard  Danger. 3 H301 Toxic if swallowed. 3 H311 Toxic in contact with skin. 3 H331 Toxic if inhaled.	≤ 1

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	May react explosively when exposed to heat. Flammable.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Methyl alcohol:
-------------------------------	--

	<p>NIOSH REL: TWA 200 ppm (260 mg/m³) ST 250 ppm (325 mg/m³) [skin] OSHA PEL: TWA 200 ppm (260 mg/m³) ACGIH TLV: TWA 200 ppm (262 mg/m³) STEL 250 ppm (328 mg/m³)</p> <p>Acetic acid, glacial; ACGIH TLV: 15ppm (37 mg/m³) STEL 10 ppm (25 mg/m³) TWA NIOSH REL: 15 ppm (37 mg/m³) STEL 10 ppm (25 mg/m³) TWA OSHA PEL: 10 ppm (25 mg/m³) TWA</p> <p>Formaldehyde: NIOSH REL: Ca TWA 0.016 ppm C 0.1 ppm OSHA PEL: [1910.1048] TWA 0.75 ppm ST 2 ppm ACGIH TLV: 0.3 ppm C</p>
8.2 Exposure controls	<p>Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
Engineering Controls	<p>Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.</p>

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Translucent yellow.
Odor	Strong formaldehyde odor.
Odor Threshold	Unknown
pH	2.0 ± 1.5
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown

Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Acids, anhydrides, aniline, bases, isocyanates, metals, oxidizing agents, phenols, reducing agents.
10.6 Hazardous Decomposition Materials	Oxides of carbon and nitrogen.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Corrosive to skin and mucous membranes. Serious Eye Damage/Irritation. Corrosive to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information


12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.

12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: 3265 ADR/RID: 3265 IMDG: 3265 IATA: 3265
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Corrosive Liquid, Acidic, Organic, N.O.S. (Formaldehyde, Acetic acid) IMDG: Corrosive Liquid, Acidic, Organic, N.O.S. (Formaldehyde, Acetic acid) IATA: Corrosive Liquid, Acidic, Organic, N.O.S. (Formaldehyde, Acetic acid) DOT: Corrosive Liquid, Acidic, Organic, N.O.S. (Formaldehyde, Acetic acid)
14.3 Transport Hazard Class(es)	DOT: 8 IATA: 8 IMDG: 8 ADR/RID: 8 
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: III IMDG: III IATA: III DOT: III
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous	None of the components in this mixture are listed.

Substances; Section 355	
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	Methanol is listed as a having developmental toxicity.
Right To Know Components	<p>Listed on Right to Know: Picric Acid, CAS# 88-89-1 California Massachusetts Minnesota New Jersey Pennsylvania Rhode Island</p> <p>Listed on Right to Know: Acetic Acid, CAS# 64-19-7 California Massachusetts Minnesota New Jersey Pennsylvania Rhode Island</p> <p>Listed on Right to Know: Formaldehyde, CAS# 50-00-0 California Massachusetts Minnesota New Jersey Pennsylvania Rhode Island</p> <p>Listed on Right to Know: Methanol, CAS# 67-56-1 California Massachusetts Minnesota New Jersey Pennsylvania Rhode Island</p>
Other Regulations	Picric Acid, CAS# 88-89-1 is required to be identified under SARA sections 302/304, 313, CERCLA, TSCA, and/or require an OSHA safety plan.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	<p>R4 EN: Forms very sensitive explosive metallic compounds – when dried. DA: Danner meget følsomme eksplosive metalforbindelser. DE: Bildet hochempfindliche explosionsgefährliche Metallverbindungen. EL: Σχηματίζει πολύ ευαίσθητες εκρηκτικές μεταλλικές ενώσεις. ES: Forma compuestos metálicos explosivos muy sensibles. FR: Forme des composés métalliques explosifs tres sensibles. IT: Forma composti metallici esplosivi molto sensibili. NL: Vormt met metalen zeer gemakkelijk ontplofbare verbindingen. PT: Forma compostos metálicos explosivos muito sensíveis. SV: Bildar mycket känsliga explosiva metallföreningar.</p> <p>R20/21/22 EN: Harmful by inhalation, in contact with skin and if swallowed. DA: Farlig ved indånding, ved hudkontakt og ved indtagelse. DE: Gesundheitsschädlich beim Einatmen, Verschlucken und Berührung mit der Haut. EL: Επιβλαβές όταν εισπνέεται, σε επαφή με το δέρμα και σε περίπτωση καταπόσεως.</p>

	<p>ES: Nocivo por inhalación, por ingestión y en contacto con la piel. FR: Nocif par inhalation, par contact avec la peau et par ingestion. IT: Nocivo per inalazione, contatto con la pelle e per ingestione. NL: Schadelijk bij inademing, opname door de mond en aanraking met de huid. PT: Nocivo por inalação, em contacto com a pele e por ingestão. SV: Farligt vid inandning, hudkontakt och förtäring.</p> <p>R45 EN: May cause cancer. DE: Kann Krebs erzeugen. FR: Peut provoquer le cancer. IT: Può provocare il cancro. ES: Puede causar cáncer. FR: Risque de lésions oculaires graves. IT: Rischio di gravi lesioni oculari. NL: Gevaar voor ernstig oogletsel. PT: Risco de lesões oculares graves. SV: Risk för allvarliga ögonskador.</p>
<p>S phrases</p>	<p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen. PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista. SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S45 EN: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). DA: Ved ulykkestilfælde eller ved ildebefindende er omgående lægebehandling nødvendig; vis etiketten, hvis det er muligt. DE: Bei Unfall oder Unwohlsein sofort Arzt zuziehen (wenn möglich, dieses Etikett vorzeigen). EL: Σε περίπτωση ατυχήματος ή αν αισθανθείτε αδιαθεσία ζητήστε αμέσως ιατρική συμβουλή (δείξτε την επικήτα αν είναι δυνατό). ES: En caso de accidente o malestar, acúdase inmediatamente al médico (si es posible, muéstresele la etiqueta). FR: En cas d'accident ou de malaise, consulter immédiatement un medecin (si possible lui montrer l'etiquette). IT: In caso di incidente o di malessere consultare immediatamente il medico (se possibile, mostrargli l'etichetta). NL: Bij een ongeval of indien men zich onwel voelt, onmiddellijk een arts raadplegen (indien mogelijk hem dit etiket tonen). PT: Em caso de acidente ou de indisposição, consultar imediatamente o médico (se possível mostrar-lhe o rótulo). SV: Vid olycksfall, illamående eller annan påverkan, kontakta omedelbart läkare. Visa om möjligt etiketten.</p> <p>S51</p>

	<p>EN: Use only in well-ventilated areas. DA: Må kun bruges på steder med god ventilation. DE: Nur in gut gelüfteten Bereichen verwenden. EL: Να χρησιμοποιείται μόνο σε καλά αεριζόμενο χώρο. ES: Útese únicamente en lugares bien ventilados. FR: Utiliser seulement dans des zones bien ventilées. IT: Usare soltanto in luogo ben ventilato. NL: Uitsluitend op goed geventileerde plaatsen gebruiken. PT: Utilizar somente em locais bem ventilados. SV: Sörj för god ventilation.</p> <p>S60 EN: This material and its container must be disposed of as hazardous waste. DA: Undgå udledning til miljøet. Se særlig vejledning/leverandørbrugsanvisning. DE: Dieses Produkt und sein Behälter sind als gefährlicher Abfall zu entsorgen. EL: Αποφύγετε την ελευθέρωσή του στο περιβάλλον. Αναφερθείτε σε ειδικές οδηγίες/ελτίο δεδομένων ασφαλείας. ES: Elimínense el producto y su recipiente como residuos peligrosos. FR: Éliminer le produit et son récipient comme un déchet dangereux. IT: Questo materiale e il suo contenitore devono essere smaltiti come rifiuti pericolosi. NL: Deze stof en de verpakking als gevaarlijk afval afvoeren. PT: Este produto e o seu recipiente devem ser eliminados como resíduos perigosos. SV: Detta material och dess behållare skall tas om hand som farligt avfall.</p>
--	--

Section 16. Other Information

Revision Date: 2019-01-11


Notice to Reader:




To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Hematoxylin, Weigert's Iron (Part A) Product Number: SSC1060	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 GHS Classification	Flammable liquids (Category 2) – H225 Skin irritation (Category 2) – H315 Eye irritation (Category 2A) – H319 Specific target organ toxicity (Category 3), Central nervous system – H336
Classification: Directive 1999/45/EC	 <p>F R11, R36; S7, S16, S24/25 S26 EN: Highly flammable DA: Meget brandfarlig DE: Leichtentzündlich EL: Πολύ εύφλεκτο ES: Fácilmente inflamable FR: Facilement inflammable IT: Facilmente infiammabile NL: Licht ontvlambaar PT: Facilmente inflamável SV: Mycket brandfarligt</p> For full text of R- & S- phrases: see Section 15.




2.2 Label Elements	Pictogram(s):   Signal word: Danger Hazard statement(s): H226 – Flammable liquid and vapour H319 – Causes serious eye irritation H315 – Causes skin irritation H336 – May cause drowsiness or dizziness Precautionary statement(s): P210 – Keep away from heat/sparks/open flames/hot surfaces – No smoking P302+ P352 – IF ON SKIN: Wash with plenty of soap and water P305 + P351 – IF IN EYES: Rinse cautiously with water for several minutes.								
NFPA Scale: 0-4									
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr> <td>HEALTH</td> <td>1</td> </tr> <tr> <td>FLAMMABILITY</td> <td>3</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>G</td> </tr> </table>	HEALTH	1	FLAMMABILITY	3	PHYSICAL HAZARD	0	PERSONAL PROTECTION	G
HEALTH	1								
FLAMMABILITY	3								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	G								
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.								


Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Ethanol	64-17-5	200-578-6	 Danger. 2 H225 Highly flammable liquid and vapour.	≤ 95
Isopropanol	67-63-0	200-661-7	 Warning. 2. H319 Causes serious eye irritation. 3. H336 May cause drowsiness or dizziness.  Danger. 3 H225 Highly flammable liquid and vapour.	≤ 5

Hematoxylin	517-28-2	208-237-3	 Warning. 4 H302 Harmful if swallowed. 2 H315 Causes skin irritation. 2 H319 Causes serious eye irritation. 3 H335 May cause respiratory irritation.	≤ 1
-------------	----------	-----------	---	-----

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water. Get medical attention if irritation continues.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosive hazards expected
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: <i>Ethanol:</i> NIOSH REL: 1000 ppm (1900 mg/m ³) TWA OSHA PEL: 1000 ppm (1900 mg/m ³) TWA ACGIH TLV: 1000ppm (1900mg/m ³) TWA <i>Isopropanol:</i> NIOSH REL: TWA 400 ppm (980 mg/m ³) ST 500 ppm (1225 mg/m ³) OSHA PEL: TWA 400 ppm (980 mg/m ³) ACGIH TLV: TWA 200 ppm STEL 400 ppm
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Transparent red/brown
Odor	Mild alcohol odor
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown

Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Strong oxidizers, silver salts, acid chlorides, alkali metals, metal hydrides, hydrazine.
10.6 Hazardous Decomposition Materials	Will not occur.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Irritating to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: 1170 ADR/RID: 1170 IMDG: 1170 IATA: 1170
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Ethanol Solutions IMDG: Ethanol Solutions IATA: Ethanol Solutions DOT: Ethanol Solutions
14.3 Transport Hazard Class(es)	DOT: Class 3 IATA: Class 3 IMDG: Class 3 ADR/RID: Class 3
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: II IMDG: II IATA: II DOT: II
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	<u>Ethanol</u> Massachusetts Pennsylvania New jersey <u>Isopropanol</u> Massachusetts Pennsylvania New jersey <u>Hematoxylin</u> Pennsylvania New Jersey
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R11

	<p>EN: Highly flammable. DA: Meget brandfarlig. DE: Leichtentzündlich. EL: Πολύ εύφλεκτο. ES: Fácilmente inflamable. FR: Facilement inflammable. IT: Facilmente infiammabile. NL: Licht ontvlambaar. PT: Facilmente inflamável. SV: Mycket brandfarligt.</p> <p>R36 EN: Irritating to eyes. DE: Reizt die Augen. FR: Irritant pour les yeux. IT: Irritante per gli occhi. ES: Irrita los ojos. FR: Irritant pour les yeux. IT: Irritante per gli occhi. NL: Irriterend voor de ogen. PT: Irritante para os olhos. SV: Irriterar ögonen.</p>
<p>S phrases</p>	<p>S7 EN: Keep container tightly closed. DA: Emballagen skal holdes tæt lukket. DE: Behälter dicht geschlossen halten. EL: Το δοχείο να διατηρείται ερμητικά κλεισμένο. ES: Manténgase el recipiente bien cerrado. FR: Conserver le récipient bien fermé. IT: Conservare il recipiente ben chiuso. NL: In goed gesloten verpakking bewaren. PT: Manter o recipiente bem fechado. SV: Förpackningen förvaras väl tillsluten.</p> <p>S16 EN: Keep away from sources of ignition - No smoking. DA: Holdes væk fra antændelseskilder - Rygning forbudt. DE: Von Zündquellen fernhalten - Nicht rauchen. EL: Μακριά από πηγές ανάφλεξης - Απαγορεύεται το κάπνισμα. ES: Conservar alejado de toda llama o fuente de chispas - No fumar. FR: Conserver a l'ecart de toute flamme ou source d'étincelles - Ne pas fumer. IT: Conservare lontano da fiamme e scintille - Non fumare. NL: Verwijderd houden van ontstekingsbronnen - Niet roken. PT: Manter afastado de qualquer chama ou fonte de ignição - Não fumar. SV: Förvaras åtskilt från antändningskällor - Rökning förbjuden</p> <p>S24/25 EN: Avoid contact with skin and eyes. DA: Undgå kontakt med huden og øjnene. DE: Berührung mit den Augen und der Haut vermeiden. EL: Αποφύγετε επαφή με το δέρμα και με τα μάτια. ES: Evítase el contacto con los ojos y la piel. FR: Éviter le contact avec la peau et les yeux. IT: Evitare il contatto con gli occhi e con la pelle. NL: Aanraking met de ogen en de huid vermijden. PT: Evitar o contacto com a pele e os olhos. SV: Undvik kontakt med huden och ögonen.</p> <p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.</p>

	<p>DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen. PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista. SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p>
--	---

Section 16. Other Information

Revision Date: 2019-01-11


Notice to Reader:



To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Hematoxylin, Weigert's Iron (Part B) Product Number: SSC1061	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik.h PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 GHS Classification	Skin irritation (Category 2) – H315 Serious eye damage (Category 1) – H318	
Classification: Directive 1999/45/EC		C R22, R41; S24/25, S28, S37/39 EN: Corrosive DA: Ætsende DE: Ätzend EL: Διαβρωτικό ES: Corrosivo FR: Corrosif IT: Corrosivo NL: Bijtend PT: Corrosivo SV: Frätande
	For full text of R- & S- phrases: see Section 15.	






2.2 Label Elements	Pictogram(s):  Signal word: Danger Hazard statement(s): H315 – Causes skin irritation H320 – Causes serious eye damage Precautionary statement(s): P302+ P352 – IF ON SKIN: Wash with plenty of soap and water P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Get immediate medical advice/attention.			
NFPA Scale: 0-4				
HMIS (U.S.A.) Scale: 0-4	HEALTH	2		
	FLAMMABILITY	0		
	PHYSICAL HAZARD	0		
	PERSONAL PROTECTION	C		
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.			

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Ferric Chloride	10025-77-1	231-729-4	 Warning. 1. H290 May be corrosive to metals.  Warning. 4 H302 Harmful if swallowed. 2 H315 Causes skin irritation.  Danger. 1 H318 Causes serious eye damage.	≤ 1
Hydrochloric Acid (~37%)	7647-01-0	231-595-7	 Warning. 3. H335 May cause respiratory irritation.  Danger. 1B. H314 Causes severe skin burns and eye damage.	≤ 1

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Ferric chloride: TWA 1 mg/m3 STEL 2 mg/m3
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Orange
Odor	Odorless
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Metals, allyl chloride, sodium, potassium.
10.6 Hazardous Decomposition Materials	Hydrogen Chloride Fumes.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Causes serious damage to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods. ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods.
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods. IATA: Not dangerous goods. IMDG: Not dangerous goods. ADR/RID: Not dangerous goods.
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R22 EN: Harmful if swallowed. DA: Farlig ved indtagelse. DE: Gesundheitsschädlich beim Verschlucken. EL: Επιβλαβές σε περίπτωση καταπόσεως. ES: Nocivo por ingestión. FR: Nocif en cas d'ingestion. IT: Nocivo per ingestione. NL: Schadelijk bij opname door de mond. PT: Nocivo por ingestão. SV: Farligt vid förtäring. R41 EN: Risk of serious damage to eyes. DA: Risiko for alvorlig øjensskade. DE: Gefahr ernster Augenschäden. EL: Κίνδυνος σοβαρών οφθαλμικών βλαβών. ES: Riesgo de lesiones oculares graves. FR: Risque de lésions oculaires graves. IT: Rischio di gravi lesioni oculari. NL: Gevaar voor ernstig oogletsel. PT: Risco de lesões oculares graves. SV: Risk för allvarliga ögonskador.

S phrases	<p>S24/25 EN: Avoid contact with skin and eyes. DA: Undgå kontakt med huden og øjnene. DE: Berührung mit den Augen und der Haut vermeiden. EL: Αποφεύγετε επαφή με το δέρμα και με τα μάτια. ES: Evítese el contacto con los ojos y la piel. FR: Éviter le contact avec la peau et les yeux. IT: Evitare il contatto con gli occhi e con la pelle. NL: Aanraking met de ogen en de huid vermijden. PT: Evitar o contacto com a pele e os olhos. SV: Undvik kontakt med huden och ögonen.</p> <p>S28 EN: After contact with skin, wash immediately with plenty of water. DA: Kommer stof på huden vaskes straks med store mængder ... (angives af fabrikanten). DE: Bei Berührung mit der Haut sofort abwaschen mit viel ... (vom Hersteller anzugeben). EL: Σε περίπτωση επαφής με το δέρμα, πλυθείτε αμέσως με άφθονο ... (το είδος του υγρού καθορίζεται από τον παραγωγό). ES: En caso de contacto con la piel, lávese inmediata y abundantemente con ... (productos a especificar por el fabricante). FR: Après contact avec la peau, se laver immédiatement et abondamment avec ... (produits appropriés à indiquer par le fabricant). IT: In caso di contatto con la pelle lavarsi immediatamente ed abbondantemente con ... (prodotti idonei da indicarsi da parte del fabbricante). NL: Na aanraking met de huid onmiddellijk wassen met veel ... (aan te geven door de fabrikant). PT: Após contacto com a pele, lavar imediata e abundantemente com ... (produtos adequados a indicar pelo produtor). SV: Vid kontakt med huden tvätta genast med mycket ... (anges av tillverkaren).</p> <p>S37/39 EN: Wear suitable gloves and eye/face protection. DA: Brug egnede beskyttelseshandsker og -briller/ansigtsskærm under arbejdet. DE: Bei der Arbeit geeignete Schutzhandschuhe und Schutzbrille/Gesichtsschutz tragen. EL: Φοράτε κατάλληλα γάντια και συσκευή προστασίας ματιών/προσώπου. ES: Úsense guantes adecuados y protección para los ojos/la cara. FR: Porter des gants appropriés et un appareil de protection des yeux/du visage. IT: Usare guanti adatti e proteggersi gli occhi/la faccia. NL: Draag geschikte handschoenen en een beschermingsmiddel voor de ogen/het gezicht. PT: Usar luvas e equipamento protector para os olhos /face adequados. SV: Använd lämpliga skyddshandskar samt skyddsglasögon eller ansiktsskydd.</p>
------------------	---

Section 16. Other Information

Revision Date: 2019-01-11


Notice to Reader:


To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Martius Yellow Solution Product Number: SSC1283	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703 U.S.A.
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com
1.4 Emergency Telephone	Chemtrec (USA): 1-800-424-9300	

Section 2. Hazards Identification

2.1 GHS Classification	Serious eye damage (Category 1) – H318 Skin irritation (Category 2) – H315 Flammable Liquids (Category 2) – H225	
2.2 Label Elements	Pictogram(s): 	
	Signal word:	Danger
	Hazard statement(s):	H303 – May be harmful if swallowed H315 – Causes skin irritation H320 – Causes serious eye damage H225 – Highly flammable liquid and vapor





	Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes. P337 + P315 – If eye irritation persists get immediate medical advice/attention. Precautionary statement(s) not listed on label: P301 + P312 - IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. <small>See Annex III and IV of the EU CLP regulation for translations of statements</small>								
NFPA Scale: 0-4 (Estimated for Mixtures and Kits)									
HMIS (U.S.A.) Scale: 0-4 (Estimated for Mixtures and Kits)	<table border="1"> <tr> <td>HEALTH</td> <td>1</td> </tr> <tr> <td>FLAMMABILITY</td> <td>3</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>H</td> </tr> </table>	HEALTH	1	FLAMMABILITY	3	PHYSICAL HAZARD	0	PERSONAL PROTECTION	H
HEALTH	1								
FLAMMABILITY	3								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	H								
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.								

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Ethanol	64-17-5	200-578-6	 Danger. 2 H225 Highly flammable liquid and vapour.	≤ 95
Phosphotungstic acid	12501-23-4	n/a	 Danger. 1B, H314 Causes severe skin burns and eye damage.	≤ 2
Isopropyl Alcohol	67-63-0	200-661-7	 Warning. 2. H319 Causes serious eye irritation. 3. H336 May cause drowsiness or dizziness.  Danger. 3 H225 Highly flammable liquid and vapour.	≤ 4

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.
Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.
Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.
4.2 Most important symptoms and effects, both acute and delayed
 See section 2.2 and 11.
4.3 Indication of any immediate medical attention and special treatment needed
 No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosive hazards expected
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures
 Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.
6.2 Environmental precautions
 Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.
6.3 Methods and materials for containment and cleaning up
 Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.
 Avoid contact with skin and eyes.
 Wash thoroughly after handling.
 Avoid breathing vapor.
7.2 Conditions for safe storage, including any incompatibilities.
 Store in well ventilated area.
 Keep container tightly closed.
 Store at 15-30°C.
7.3 Specific end use(s).
 See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Ethanol: NIOSH REL: 1000 ppm (1900 mg/m ³) TWA OSHA PEL: 1000 ppm (1900 mg/m ³) TWA ACGIH TLV: 1000ppm (1900mg/m ³) TWA
-------------------------------	--

	<p>Isopropyl alcohol: NIOSH REL: TWA 400 ppm (980 mg/m³) ST 500 ppm (1225 mg/m³) OSHA PEL: TWA 400 ppm (980 mg/m³) ACGIH TLV: TWA 200 ppm STEL 400 ppm</p>
8.2 Exposure controls	<p>Personal Protective Equipment (PPE):</p> <p>Eye/Face protection. Safety glasses or goggles are required.</p> <p>Skin protection. Protective clothing is required.</p> <p>Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove.</p> <p>Respiratory protection. Avoid breathing vapor.</p> <p>Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
Engineering Controls	<p>Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.</p>

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Yellow
Odor	Mild alcohol odor.
Odor Threshold	Unknown
pH	Acidic
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity



10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Metals, strong bases, acetyl chloride, acetyl bromide, nitric acid, iodine-mercuric oxide, hydrogen-palladium, chromic anhydride, magnesium.
10.6 Hazardous Decomposition Materials	Hydrogen chloride, oxides of carbon.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	<p>Acute Toxicity. No relevant data available</p> <p>Skin Corrosion/Irritation. Irritating to skin and mucous membranes.</p> <p>Serious Eye Damage/Irritation. Seriously irritating to eye.</p> <p>Respiratory or skin sensitization. No relevant data available.</p> <p>Germ Cell Mutagenicity. No relevant data available.</p>
Carcinogenicity.	<p>International Agency for Research on Cancer (IARC). None of the components are listed.</p> <p>National Toxicology Program (NTP). None of the components are listed.</p>


Section 12. Ecological Information

12.1 Toxicity	<p>Fish: No relevant studies identified.</p> <p>Crustacea: No relevant studies identified.</p> <p>Algae/Aquatic Plants: No relevant studies identified.</p> <p>Other Organisms: No relevant studies identified.</p>
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	<p>PBT: This mixture does not contain any substances that are assessed to be a PBT.</p> <p>vPvB: This mixture does not contain any substances that are assessed to be a vPvB.</p>

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information – DOT IATA, IMDG, ADR/RID

14.1 UN Number	1170
14.2 UN Proper Shipping Name	Ethyl alcohol solution
14.3 Transport Hazard Class(es)	Class 3 
14.4 Packing Group DOT, IATA, IMDG, ADR	II
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	<i>Ethanol</i> Massachusetts New Jersey Pennsylvania <i>Isopropyl Alcohol</i> Massachusetts New Jersey Pennsylvania


Section 16. Other Information


The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Cancer Diagnostics, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

Section 1. Identification

1.1 Product Identifier	Product Name: Crystal Scarlet Solution Product Number: SSC1284	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703 U.S.A.
	Phone Number	877-846-5393
	Fax Number	877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com
1.4 Emergency Telephone	Chemtrec (USA): 1-800-424-9300	

Section 2. Hazards Identification

2.1 GHS Classification	Skin irritation (Category 2) – H315 Eye irritation (Category 2) – H320	
2.2 Label Elements	Pictogram(s):	
	Signal word:	Warning
	Hazard statement(s):	H315 – Causes skin irritation H320 – Causes eye irritation
	Precautionary statement(s):	P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes.



NFPA Scale: 0-4 (Estimated for Mixtures and Kits)			
HMIS (U.S.A.) Scale: 0-4 (Estimated for Mixtures and Kits)	HEALTH	1	
	FLAMMABILITY	0	
	PHYSICAL HAZARD	0	
	PERSONAL PROTECTION	B	
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.		

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Acetic Acid	64-19-7	200-580-7	 Warning. 3, H226 Flammable liquid and vapour.  Danger. 1A H314 Causes severe skin burns and eye damage.	≤ 1

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

<p>8.1 Control parameters</p>	<p>Exposure Limits: Acetic acid, glacial; ACGIH TLV: 15ppm (37 mg/m³) STEL 10 ppm (25 mg/m³) TWA NIOSH REL: 15 ppm (37 mg/m³) STEL 10 ppm (25 mg/m³) TWA OSHA PEL: 10 ppm (25 mg/m³) TWA</p>
<p>8.2 Exposure controls</p>	<p>Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
<p>Engineering Controls</p>	<p>Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.</p>



Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Dark Red
Odor	Slight vinegar smell.
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Oxidizing agents, Strong oxidizing agents, Strong acids, Metals, Amines, Alcohols, Peroxides, permanganates, Soluble carbonates and phosphates, Hydroxides.
10.6 Hazardous Decomposition Materials	Carbon oxides.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	<p>Acute Toxicity. No relevant data available</p> <p>Skin Corrosion/Irritation. Irritating to skin and mucous membranes.</p> <p>Serious Eye Damage/Irritation. Irritating to eye.</p> <p>Respiratory or skin sensitization.</p>
--	---

	No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number	Not regulated for transport.
14.2 UN Proper Shipping Name	Not regulated for transport.
14.3 Transport Hazard Class(es)	Not regulated for transport.
14.4 Packing Group	Not regulated for transport.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.



Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.


Section 16. Other Information


The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Cancer Diagnostics, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

Section 1. Identification

1.1 Product Identifier	Product Name: Phosphotungstic Acid Differentiator Product Number: SSC1285	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancerdiagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703 U.S.A.
	Phone Number	877-846-5393
	Fax Number	788-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com
1.4 Emergency Telephone	Chemtrec (USA): 1-800-424-9300	

Section 2. Hazards Identification

2.1 GHS Classification	Serious eye damage (Category 1) – H318 Skin irritation (Category 2) – H315 Flammable Liquids (Category 2) – H225	
2.2 Label Elements	Pictogram(s): 	
	Signal word:	Danger
	Hazard statement(s):	H303 – May be harmful if swallowed H315 – Causes skin irritation H320 – Causes serious eye damage H225 – Highly flammable liquid and vapor





	Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes. P337 + P315 – If eye irritation persists get immediate medical advice/attention. Precautionary statement(s) not listed on label: P301 + P312 - IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. <small>See Annex III and IV of the EU CLP regulation for translations of statements</small>								
NFPA Scale: 0-4 (Estimated for Mixtures and Kits)									
HMIS (U.S.A.) Scale: 0-4 (Estimated for Mixtures and Kits)	<table border="1"> <tr> <td>HEALTH</td> <td>1</td> </tr> <tr> <td>FLAMMABILITY</td> <td>3</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>H</td> </tr> </table>	HEALTH	1	FLAMMABILITY	3	PHYSICAL HAZARD	0	PERSONAL PROTECTION	H
HEALTH	1								
FLAMMABILITY	3								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	H								
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.								

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%/w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Ethanol	64-17-5	200-578-6	 Danger. 2 H225 Highly flammable liquid and vapour.	≤ 60
Phosphotungstic acid	12501-23-4	n/a	 Danger. 1B, H314 Causes severe skin burns and eye damage.	≤ 2
Isopropyl Alcohol	67-63-0	200-661-7	 Warning. 2. H319 Causes serious eye irritation. 3. H336 May cause drowsiness or dizziness.  Danger. 3 H225 Highly flammable liquid and vapour.	≤ 4

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.
Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.
Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.
4.2 Most important symptoms and effects, both acute and delayed
 See section 2.2 and 11.
4.3 Indication of any immediate medical attention and special treatment needed
 No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosive hazards expected
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures
 Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.
6.2 Environmental precautions
 Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.
6.3 Methods and materials for containment and cleaning up
 Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.
 Avoid contact with skin and eyes.
 Wash thoroughly after handling.
 Avoid breathing vapor.
7.2 Conditions for safe storage, including any incompatibilities.
 Store in well ventilated area.
 Keep container tightly closed.
 Store at 15-30°C.
7.3 Specific end use(s).
 See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Ethanol: NIOSH REL: 1000 ppm (1900 mg/m ³) TWA OSHA PEL: 1000 ppm (1900 mg/m ³) TWA ACGIH TLV: 1000ppm (1900mg/m ³) TWA
-------------------------------	--



	<p>Isopropyl alcohol: NIOSH REL: TWA 400 ppm (980 mg/m³) ST 500 ppm (1225 mg/m³) OSHA PEL: TWA 400 ppm (980 mg/m³) ACGIH TLV: TWA 200 ppm STEL 400 ppm</p>
8.2 Exposure controls	<p>Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
Engineering Controls	<p>Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.</p>

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Colorless
Odor	Mild alcohol odor.
Odor Threshold	Unknown
pH	Acidic
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Metals, strong bases, acetyl chloride, acetyl bromide, nitric acid, iodine-mercuric oxide, hydrogen-palladium, chromic anhydride, magnesium.
10.6 Hazardous Decomposition Materials	Hydrogen chloride, oxides of carbon.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	<p>Acute Toxicity. No relevant data available</p> <p>Skin Corrosion/Irritation. Irritating to skin and mucous membranes.</p> <p>Serious Eye Damage/Irritation. Seriously irritating to eye.</p> <p>Respiratory or skin sensitization. No relevant data available.</p> <p>Germ Cell Mutagenicity. No relevant data available.</p>
Carcinogenicity.	<p>International Agency for Research on Cancer (IARC). None of the components are listed.</p> <p>National Toxicology Program (NTP). None of the components are listed.</p>


Section 12. Ecological Information

12.1 Toxicity	<p>Fish: No relevant studies identified.</p> <p>Crustacea: No relevant studies identified.</p> <p>Algae/Aquatic Plants: No relevant studies identified.</p> <p>Other Organisms: No relevant studies identified.</p>
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	<p>PBT: This mixture does not contain any substances that are assessed to be a PBT.</p> <p>vPvB: This mixture does not contain any substances that are assessed to be a vPvB.</p>

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information – DOT IATA, IMDG, ADR/RID

14.1 UN Number	1170
14.2 UN Proper Shipping Name	Ethyl alcohol solution
14.3 Transport Hazard Class(es)	Class 3 
14.4 Packing Group DOT, IATA, IMDG, ADR	II
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	<i>Ethanol</i> Massachusetts New Jersey Pennsylvania <i>Isopropyl Alcohol</i> Massachusetts New Jersey Pennsylvania

Section 16. Other Information


Revision Date: 2022-12-22



The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Cancer Diagnostics, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product.

Section 1. Identification

1.1 Product Identifier	Product Name: Aniline Blue (MSB) Product Number: SSC1285	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 GHS Classification	Skin irritation (Category 2) – H315 Eye irritation (Category 2A) – H319	
Classification: Directive 1999/45/EC	 Xi, R36/38; S24/25, S37/39 EN: Irritant DA: Lokalirriterende DE: Reizend EL: Ερεθιστικό ES: Irritante FR: Irritant IT: Irritante NL: Irriterend PT: Irritante SV: Irriterande For full text of R- & S- phrases: see Section 15.	



2.2 Label Elements	Pictogram(s):  Signal word: Warning Hazard statement(s): H315 – Causes skin irritation H319 – Causes serious eye irritation Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes.		
NFPA Scale: 0-4			
HMIS (U.S.A.) Scale: 0-4	HEALTH	1	
	FLAMMABILITY	0	
	PHYSICAL HAZARD	0	
	PERSONAL PROTECTION	B	
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.		


Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Acetic Acid	64-19-7	200-580-7	 Warning. 3, H226 Flammable liquid and vapour.  Danger. 1A H314 Causes severe skin burns and eye damage.	≤ 2

Aniline Blue	28631-66-5	249-113-9	 <p>Warning. 2 H315 Causes skin irritation. 1, H317 May cause an allergic skin reaction. 2 H319 Causes serious eye irritation</p>	≤ 2
--------------	------------	-----------	--	-----

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water. Get immediate medical attention if irritation persists.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	Carbon oxides, Nitrogen oxides (NO _x), Sulphur oxides, Sodium oxides.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage



<p>7.1 Precautions for safe handling. Avoid contact with skin and eyes. Wash thoroughly after handling. Avoid breathing vapor.</p> <p>7.2 Conditions for safe storage, including any incompatibilities. Store in well ventilated area. Keep container tightly closed. Store at 15-30°C.</p> <p>7.3 Specific end use(s). See section 1.2</p>

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	<p>Exposure Limits: Acetic acid, glacial; ACGIH TLV: 15ppm (37 mg/m³) STEL 10 ppm (25 mg/m³) TWA NIOSH REL: 15 ppm (37 mg/m³) STEL 10 ppm (25 mg/m³) TWA OSHA PEL: 10 ppm (25 mg/m³) TWA</p>
8.2 Exposure controls	<p>Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
Engineering Controls	<p>Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.</p>

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Dark Blue
Odor	Odorless
Odor Threshold	Unknown

pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Strong oxidizing agents, Bases, Acids, Metals.
10.6 Hazardous Decomposition Materials	Carbon monoxide, carbon dioxide, nitrogen oxides.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Seriously irritating to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP).

	None of the components are listed.
--	------------------------------------

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods. ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods.
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods. IATA: Not dangerous goods. IMDG: Not dangerous goods. ADR/RID: Not dangerous goods.
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	<u>Acetic Acid, CAS # 64-19-7</u> California Massachusetts Minnesota New Jersey Pennsylvania Rhode Island
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R36/38 EN: Irritating to eyes and skin. DA: Irriterer øjnene og huden DE: Reizt die Augen und die Haut. EL: Ερεθίζει τα μάτια και το δέρμα. ES: Irrita los ojos y la piel. FR: Irritant pour les yeux et la peau. IT: Irritante per gli occhi e la pelle. NL: Irriterend voor de ogen en de huid. PT: Irritante para os olhos e pele. SV: Irriterar ögonen och huden.
S phrases	S24/25 EN: Avoid contact with skin and eyes. DA: Undgå kontakt med huden og øjnene. DE: Berührung mit den Augen und der Haut vermeiden. EL: Αποφεύγετε επαφή με το δέρμα και με τα μάτια. ES: Evítese el contacto con los ojos y la piel. FR: Éviter le contact avec la peau et les yeux. IT: Evitare il contatto con gli occhi e con la pelle. NL: Aanraking met de ogen en de huid vermijden. PT: Evitar o contacto com a pele e os olhos. SV: Undvik kontakt med huden och ögonen. S37/39 EN: Wear suitable gloves and eye/face protection. DA: Brug egnede beskyttelseshandsker og -briller/ansigtsskærm under arbejdet. DE: Bei der Arbeit geeignete Schutzhandschuhe und Schutzbrille/ Gesichtsschutz tragen. EL: Φοράτε κατάλληλα γάντια και συσκευή προστασίας ματιών/προσώπου. ES: Úsense guantes adecuados y protección para los ojos/la cara. FR: Porter des gants appropriés et un appareil de protection des yeux/du visage. IT: Usare guanti adatti e proteggersi gli occhi/la faccia. NL: Draag geschikte handschoenen en een beschermingsmiddel voor de ogen/het gezicht. PT: Usar luvas e equipamento protector para os olhos /face adequados. SV: Använd lämpliga skyddshandskar samt skyddsglasögon eller ansiktsskydd.



Safety Data Sheet

Emergency 1-800-424-9300
Contact (intl.) 1-352-323-3500

--	--

Section 16. Other Information

Revision Date: 2022-12-22

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.