



STAINLESS



Stainless Steel Base Molds

Stainless Steel Base Molds

Durable reusable heavy duty stainless steel base molds designed to be used with all standard cassettes.

Interior shape aids in better release and easier sectioning of the formed wax block.

COATED STAINLESS



Atomic™ Coated Stainless Steel Base Molds

ARCHIMEDES™ Series Atomic™ Coated Stainless Steel Base Mold

Archimedes™ Atomic™ version adds a permanent non-stick coating for faster cooling, seamless mold release and easier cleaning.

Stainless Steel Base Molds	Size	Qty
<input type="checkbox"/> MBM0707	7 x 7 x 5 mm	BX/12
<input type="checkbox"/> MBM1515	15 x 15 x 5 mm	BX/12
<input type="checkbox"/> MBM2424	24 x 24 x 5 mm	BX/12
<input type="checkbox"/> MBM3024	30 x 24 x 5 mm	BX/12
<input type="checkbox"/> MBM3724	37 x 24 x 5 mm	BX/12

Double-Deep Mold	Size	Qty
<input type="checkbox"/> MBM3724-XL	37 x 24 x 10 mm	BX/12

Archimedes™ Atomic™ SS Base Molds	Size	Qty
<input type="checkbox"/> MBM0707-ATM	7 x 7 x 5 mm	BX/12
<input type="checkbox"/> MBM1515-ATM	15 x 15 x 5 mm	BX/12
<input type="checkbox"/> MBM2424-ATM	24 x 24 x 5 mm	BX/12
<input type="checkbox"/> MBM3024-ATM	30 x 24 x 5 mm	BX/12
<input type="checkbox"/> MBM3724-ATM	37 x 24 x 5 mm	BX/12

For InVitro Diagnostics Use



4300 Emperor Blvd - 400 Durham, NC 27703
 Tel: 1-877-846-5393 Fax: 1-877-817-1716
cancerdiagnostics.com