
Cotton Tipped Applicator

Double tapered mini-tip, Wood handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Tapered mini cotton tip ideal for hard-to-reach or small areas. Available non-sterile or sterile.

Dimensions

5.875" Overall Length (149.23 mm)

Tip

- 0.120" Dia. (3.05 mm) Tapered to 0.050" (1.27mm)
- 0.5" Length (12.7 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.812" Length (147.62 mm)

Qualities

- · Medical grade quality
- · Tapered mini-tip
- · Rigid wood handle
- Made in USA

Applications

Critical Environments / Industrial / Medical

- Multi-purpose use
- · Wound care
- · General cleaning swab
- Removing epoxy residues from connectors
- Cleaning tubes, scopes, sensors, and other components

Ref #	Packaging
826 WC Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-826 5WC Sterile	10/20/5 (10 boxes per case, 20 packages per box, 5 swabs per bag)



Tel 800 321-2313 (US & Canada Only) or 207-876-3311

Web puritanmedproducts.com Email sales@puritanmedproducts.com