

**COMPONENTS****CAS#**

Deionized Water

7732-18-5

Formic Acid

64-18-8

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300*For Invitro Diagnostic Use Only***Formic Acid Bone Decalcifier****IMMUNO-4™****DIRECTIONS FOR USE:**

Immuno-4™ is a very mild formic acid decalcifier. It is formulated from semiconductor grade formic acid and highly purified reagent grade water. Immuno-4 should be used when you are decalcifying bone for immunohistochemistry. Immuno-4 preserves the antigenicity of both decalcified needle biopsies of bone marrow and large specimens for immunochemical studies.

**DANGER:** Contains dilute acid. May be corrosive to metallic parts, skin, eyes and clothing. Avoid prolonged contact. Causes severe skin burns and eye damage. Do not breathe dust/fume/gas/mist/vapors/spray. Wash face, hands and any exposed skin thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection.

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. Immediately call a POISON CENTER or doctor/physician.

**IF ON SKIN (OR HAIR):** Take off immediately all contaminated clothing. Rinse skin with water/shower. Take off all contaminated clothing and wash it before reuse.

**IF INHALED:** Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician.

**IF SWALLOWED:** Rinse mouth. DO NOT induce vomiting. Immediately call a POISON CENTER or doctor/physician.

**IN CASE OF FIRE:** Use CO<sub>2</sub>, dry chemical, or foam for extinction.

[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)

Cancer Diagnostics, Inc. 4300 Emperor Blvd., 400 Durham, NC 27703 (877) 846 5393 Fax: (877) 817 1716

REF M4038 | Gal (3.8L)

Store locked up at room temperature.  
Dispose of contents/container to an approved waste disposal plant.

59°F (15°C) 86°F (30°C)



MX1328 IMMUNO-4 | Gal AW-25867 Rev.B 2020-06-18

MADE IN USA