



Pureview® NBF

10% NBF OPTIMIZED FOR AUTOMATION

REF #FX1048 1 GAL (3.8L)

LOT



For In Vitro Diagnostics Use

59°F (15°C) — 86°F (30°C)

Date Received: _____ Date Opened: _____

PureView® NBF produces no precipitates in the presence of alcohol to let you skip the routine warm water processing flush. PureView® NBF is specifically formulated to inhibit formation of salts within tissue processor lines of all processors and can replace regular 10% NBF without any adjustments to procedures or protocols. Compatible with CDI's N-15 formalin neutralizer and all industry neutralizers.



DANGER
CONTAINS FORMALDEHYDE
HAZARD STATEMENTS

Harmful if inhaled. Causes skin irritation and serious eye damage. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction. May cause genetic defects. May cause cancer. May cause damage to organs (lungs, nose) through single, prolonged, or repeated exposure.

PRECAUTIONARY STATEMENTS - PREVENTION

Use only outdoors or in a well-ventilated area. Wash body thoroughly after handling. Wear protective gloves and clothing. Wear eye and face protection. Wear NIOSH approved respiratory protection. Contaminated work clothing must not be allowed out of the workplace. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Do not breathe dust, vapors. Do not eat, drink, or smoke when using this product.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. Immediately call a doctor.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

IF INHALED: If breathing is difficult, remove person to fresh air and keep comfortable for breathing. If experiencing respiratory symptoms: call a doctor. If exposed or concerned: get medical attention/call a doctor. Call a doctor if you feel unwell.

STORAGE & DISPOSAL

Store locked up. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703
Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

MX1336 Pureview NBF 1 Gal AWW-546 Rev.B 2020-07-01