




Acid Alcohol, 1%

REF #ST0016 | GAL (3.8L)

LOT



For In Vitro Diagnostics Use



Date Received: _____ Date Opened: _____

Formulated as a differentiator of Hematoxylin in regressive H&E staining procedures. This differentiator provides ideal differentiation of nuclear structures. Non-nuclear structures show clear staining by eosin with precise color variations by the eosin. Residual mucin staining by the hematoxylin is minimal to slight.



DANGER **HAZARD STATEMENTS**

Causes skin irritation. Causes serious eye irritation. Causes damage to organs. Highly flammable liquid and vapor.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash face, hands and any exposed skin thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection. Do not breathe dust/fume/gas/mist/vapors/spray. Do not eat, drink or smoke when using this product. Keep away from heat/sparks/open flames/hot surfaces. — No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical/ventilating/lighting/equipment. Use only non-sparking tools. Take precautionary measures against static discharge.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

IF ON SKIN: Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. If skin irritation occurs: Get medical advice/attention. Wash contaminated clothing before reuse.

IF exposed: Call a POISON CENTER or doctor/physician.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store locked up, in a well ventilated place at room temperature. Keep cool.

Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

MX1323 Acid Alcohol 1.0 1 Gal. AW-25812 Rev.B 2020-06-17