



Dermatoxylin™


Skin Hematoxylin

REF #ST0014 500mL

LOT



For In Vitro Diagnostics Use



Date Received: _____ Date Opened: _____

Dermatoxylin™ is a regressive skin hematoxylin optimized for dermatology specimens and is prepared to produce optimum staining in rapid time. Dermatoxylin™ is formulated for routine H&E Staining in the Dermatopathology Laboratory and is designed to work with Dermosin™ skin eosin, which has the appropriate amount of acid for proper counterstaining of epidermis and dermis.



DANGER HAZARD STATEMENTS

Causes skin and serious eye irritation. May cause damage to organs(kidney) through prolonged or repeated exposure. Suspected of damaging fertility of the unborn child.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves and clothing. Wear eye and face protection. Do not breathe dust, vapors. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood.

PRECAUTIONARY STATEMENTS - RESPONSE

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on this label). If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse

IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing If eye irritation persists: get medical attention. Call a doctor if you feel unwell. If exposed or concerned: get medical attention.

STORAGE & DISPOSAL

Store locked up at room temperature. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

MX1340 Dermatoxylin 500mL AW-1045 Rev.B 2020-07-02