



Harris Hematoxylin

Optimized, Mercury Free (Filtered)

REF #SH5777 1 Gallon (3.8L)

LOT



For In Vitro Diagnostics Use



Date Received: _____ Date Opened: _____

Harris Hematoxylin produces a deep, dark blue-purple stain. It is generally used in regressive methods, but may also be used progressively both in histology and cytology preparations.



DANGER HAZARD STATEMENTS

Causes skin and serious eye irritation. May cause damage to organs(kidney) through prolonged or repeated exposure. Suspected of damaging fertility of the unborn child.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves and clothing. Wear eye and face protection. Do not breathe dust, vapors. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue Rinsing. If eye irritation persists: get medical attention. Call a doctor if you feel unwell. If exposed or concerned: get medical attention.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store locked up at room temperature. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300

4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

MX1359 Harris Hematoxylin 1 Gallon AW-25816 Rev.C 2020-07-06

