



# Restore 4™ Concentrate

## Enzymatic Instrument Cleaner

To make working solution;  
Add 500mL to 3285mL of water

REF #REI004 | 1 Gallon (3.8L) CS/4

LOT



For In Vitro Diagnostics Use



Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_

Restore 4™ is an enzymatic detergent instrument cleaner. 4 enzymes efficiently and systematically break down proteinaceous bioburden. The detergent combination provides a light foaming action to aid in cleaning. Enzymes include; Lipase (for fat), Amylase (for starch), Carbohydrase (for high starches), and Protease (for blood)



### **DANGER** **HAZARD STATEMENTS**

Causes skin and eye irritation. Corrosive to skin. Toxic to aquatic life with long lasting effects. Causes serious eye damage

### **PRECAUTIONARY STATEMENTS - PREVENTION**

Do not breathe dust/fume/gas/mist/vapors/spray. Wash face, hands and any exposed skin thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection

### **PRECAUTIONARY STATEMENTS - RESPONSE**

Immediately call a POISON CENTER or doctor/physician. Specific treatment (see Section 4 on Safety Data Sheet).

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician

**IF ON SKIN:** Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse.

**IF INHALED:** Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician.

**IF SWALLOWED:** Rinse mouth. DO NOT induce vomiting.

**IN CASE OF FIRE:** Use CO<sub>2</sub>, dry chemical, or foam for extinction.

### **STORAGE & DISPOSAL**

Store locked up at room temperature. Dispose of contents/container to an approved waste disposal plant.

**NOTICE:** Read Safety Data Sheet for additional information.

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

**Made in USA**

MX1610 Restore 4 Concentrate Gallons, 4 AW-5271605 Rev.C 2021-01-14