



Restore 4™ Concentrate

Enzymatic Instrument Cleaner

To make working solution;
Add 500mL to 3285mL of water

REF #REI000 1 Gallon (3.8L)

LOT



For In Vitro Diagnostics Use



Date Received: _____ Date Opened: _____

Restore 4™ is an enzymatic detergent instrument cleaner. 4 enzymes efficiently and systematically break down proteinaceous bioburden. The detergent combination provides a light foaming action to aid in cleaning. Enzymes include; Lipase (for fat), Amylase (for starch), Carbohydrase (for high starches), and Protease (for blood)



DANGER **HAZARD STATEMENTS**

Causes skin and eye irritation. Corrosive to skin. Toxic to aquatic life with long lasting effects. Causes serious eye damage

PRECAUTIONARY STATEMENTS - PREVENTION

Do not breathe dust/fume/gas/mist/vapors/spray. Wash face, hands and any exposed skin thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection

PRECAUTIONARY STATEMENTS - RESPONSE

Immediately call a POISON CENTER or doctor/physician. Specific treatment (see Section 4 on Safety Data Sheet).

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician

IF ON SKIN: Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse.

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician.

IF SWALLOWED: Rinse mouth. DO NOT induce vomiting.

IN CASE OF FIRE: Use CO₂, dry chemical, or foam for extinction.

STORAGE & DISPOSAL

Store locked up at room temperature. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

MX1609 Restore 4 Concentrate 1 Gallon AW-5271606 Rev.C 2021-01-14