



# Atomic™, Coated, Non-Stick and Traditional Series Stainless Metal Base Molds



Atomic, Non-Stick Blue Coated, Reusable Heavy Duty



Atomic, Non-Stick White Coated, Reusable Heavy Duty



Stainless Steel Base Molds, Reusable Heavy Duty Standard Size



## FEATURES (ATOMIC SERIES)

Permanent metal (308 stainless steel) base molds feature non-stick, nano-layer coating for faster cooling and clean, easy block release. Same shape as traditional stainless molds.

- Non-stick (Pearlex) Coating for clean, effortless block release and no residue
- 50% faster cooling speeds up embedding process
- Blue coating offers optimal color contrast

## ORDERING INFORMATION

MBM0707-ATMB	Atomic, Non-Stick Blue Coated, Reusable Heavy Duty	7 x 7 x 5 mm
MBM1515-ATMB	Atomic, Non-Stick Blue Coated, Reusable Heavy Duty	15 x 15 x 5 mm
MBM2424-ATMB	Atomic, Non-Stick Blue Coated, Reusable Heavy Dut	24 x 24 x 5 mm
MBM3024-ATMB	Atomic, Non-Stick Blue Coated, Reusable Heavy Duty	30 x 24 x 5 mm
MBM3724-ATMB	Atomic, Non-Stick Blue Coated, Reusable Heavy Duty	37 x 24 x 5 mm



MBM0707-ATM	Atomic, Non Stick White Coated, Reusable Heavy Duty	7 x 7 x 5 mm
MBM1515-ATM	Atomic, Non Stick White Coated, Reusable Heavy Duty	15 x 15 x 5 mm
MBM2424-ATM	Atomic, Non Stick White Coated, Reusable Heavy Duty	24 x 24 x 5 mm
MBM3024-ATM	Atomic, Non Stick White Coated, Reusable Heavy Duty	30 x 24 x 5 mm
MBM3724-ATM	Atomic, Non Stick White Coated, Reusable Heavy Duty	37 x 24 x 5 mm



MBM0707	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	7 x 7 x 5 mm
MBM1515	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	15 x 15 x 5 mm
MBM2424	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	24 x 24 x 5 mm
MBM3024	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	30 x 24 x 5 mm
MBM3724	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	37 x 24 x 5 mm



MBM3724-XL	Double Deep Base Molds, Reusable Heavy Duty Standard Size	37mm x 24mm x 10mm
------------	---	--------------------



**CANCER**  
DIAGNOSTICS, INC.

### Contact Information

Visit us online

4300 Emperor Blvd-400  
Durham, North Carolina 27703

**Toll Free:** (877) 846 5393  
**Fax:** (877) 817 1716

[cancerdiagnostics.com](http://cancerdiagnostics.com)