



# Atomic™, Coated, Non-Stick and Traditional Series Stainless Metal Base Molds



**Atomic™ Blue, Non-Stick Coated, Reusable Heavy Duty**



**Atomic™ White, Non-Stick Coated, Reusable Heavy Duty**



**Stainless Steel, Reusable Heavy Duty Standard Size**



## FEATURES

Durable reusable 308 Stainless base molds. Interior mold design helps optimize sectioning. Standard 5mm and Double 10mm depth.

Atomic™ Series feature additional non-stick, nano-layer coating for faster cooling and clean, easy block release.

- Non-Stick (Pearlex ✨) Coating for clean, effortless block release and no residue.
- 50% faster cooling speeds up embedding process.
- Atomic White or Atomic Blue offers optimal color contrast.

## ORDERING INFORMATION

MBM0707-ATMB	<b>Atomic™</b> Blue, Non-Stick Coated, Reusable Heavy Duty	7 x 7 x 5 mm
MBM1111-ATMB	<b>Atomic™</b> Blue, Non-Stick Coated, Reusable Heavy Duty	11 x 11 x 5 mm
MBM1515-ATMB	<b>Atomic™</b> Blue, Non-Stick Coated, Reusable Heavy Duty	15 x 15 x 5 mm
MBM2424-ATMB	<b>Atomic™</b> Blue, Non-Stick Coated, Reusable Heavy Duty	24 x 24 x 5mm
MBM3024-ATMB	<b>Atomic™</b> Blue, Non-Stick Coated, Reusable Heavy Duty	30 x 24 x 5 mm
MBM3724-ATMB	<b>Atomic™</b> Blue, Non-Stick Coated, Reusable Heavy Duty	37 x 24 x 5 mm



MBM0707-ATM	<b>Atomic™</b> White, Non-Stick Coated, Reusable Heavy Duty	7 x 7 x 5 mm
MBM1111-ATM	<b>Atomic™</b> White, Non-Stick Coated, Reusable Heavy Duty	11 x 11 x 5 mm
MBM1515-ATM	<b>Atomic™</b> White, Non-Stick Coated, Reusable Heavy Duty	15 x 15 x 5 mm
MBM2424-ATM	<b>Atomic™</b> White, Non-Stick Coated, Reusable Heavy Duty	24 x 24 x 5 mm
MBM3024-ATM	<b>Atomic™</b> White, Non-Stick Coated, Reusable Heavy Duty	30 x 24 x 5 mm
MBM3724-ATM	<b>Atomic™</b> White, Non-Stick Coated, Reusable Heavy Duty	37 x 24 x 5 mm



MBM0707	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	7 x 7 x 5 mm
MBM1111	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	11 x 11 x 5 mm
MBM1515	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	15 x 15 x 5 mm
MBM2424	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	24 x 24 x 5 mm
MBM3024	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	30 x 24 x 5 mm
MBM3724	Stainless Steel Base Molds, Reusable Heavy Duty Standard Size	37 x 24 x 5 mm



MBM3724-XL	Stainless Steel, Reusable Heavy Duty Double Deep Size	37 x 24 x 10mm
MBM3724-ATMB-XL	<b>Atomic™</b> Blue, Reusable Heavy Duty Double Deep Size	37 x 24 x 10mm



**CANCER**  
DIAGNOSTICS, INC.

Blue Color metal base molds are trademark of Cancer Diagnostics, Inc.

**Contact Information**

**Visit us Online**

116 Page Point Circle  
Durham, North Carolina 27703

**Toll Free:** (877) 846 5393  
**Fax:** (877) 817 1716

[cancerdiagnostics.com](http://cancerdiagnostics.com)