



# Harris Hematoxylin, Acidified

Optimized, Mercury Free (Filtered)

REF #SH5775 | 1 Gallon (3.8L)

LOT



For In Vitro Diagnostics Use



Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_

Harris Hematoxylin produces a deep, dark blue-purple stain. It is generally used in regressive methods, but may also be used progressively both in histology and cytology preparations.



## WARNING HAZARD STATEMENTS

Harmful if inhaled. Causes skin irritation. Causes serious eye irritation.

### PRECAUTIONARY STATEMENTS - PREVENTION

Avoid breathing dust/fume/gas/mist/vapors/spray. Use only outdoors or in a well-ventilated area. Wash face, hands and any exposed skin thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection.

### PRECAUTIONARY STATEMENTS - RESPONSE

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

**IF ON SKIN:** Wash with plenty of soap and water. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse.

**IF ON SKIN:** Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or doctor/physician if you feel unwell.

**IN CASE OF FIRE:** Use CO<sub>2</sub>, dry chemical, or foam for extinction.

### STORAGE & DISPOSAL

Store at room temperature. Dispose of contents/container to an approved waste disposal plant.

**NOTICE:** Read Safety Data Sheet for additional information.

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)

**Made in USA**

MX1446 Harris Hematoxylin Acidified | 1 Gallon AW-26030 Rev.B 2020-08-10