



# CytoStain™ Bluing

## Cytology Bluing Solution

REF #FX2106 1000mL

LOT



**For In Vitro Diagnostics Use**



Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_

A cytology bluing solution formulated to optimally blue both regressive and progressive papanicolaou and non gyn specimens. Gentle alkaline reagents minimize cell loss from glass.

### **NON HAZARDOUS CHEMICAL**

**This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)**

Although this product does not normally require specific hazard precautions, chemical users should always take care to minimize personnel exposure and workplace contamination.

### **STORAGE & DISPOSAL -**

Store tightly closed at room temperature. Dispose of contents/container to an approved waste disposal plant.

**NOTICE:** Read Safety Data Sheet for additional information.

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703  
Tel: 1-877-846-5393 Fax: 1-877-817-1716  
[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com) **Made in USA**

MX1819 CytoStain Bluing 1000mL AWW-5417085 Rev.C 2021-06-23