



Nail-RX™ RAPID

KERATIN SOFTENING AGENT

REF #NLR05 5 Gallon CUBE

LOT



For In Vitro Diagnostics Use



Date Received: _____ Date Opened: _____

Use pre-processing to soften nail clippings. Nail-RX™ nail prep solution softens brittle nails, making them easier to cut and process, while helping to keep the specimen from washing off the slide. Soak for 1 hour and process as usual. (Longer depending on size of specimen.)



DANGER HAZARD STATEMENTS

Harmful if swallowed. Causes severe skin burns and eye damage. May be corrosive to metals.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash face, hands and any exposed skin thoroughly after handling. Do not eat, drink or smoke when using this product. Do not breathe dust/fume/gas/mist/vapors/spray. Wear protective gloves/protective clothing/eye protection/face protection. Keep only in original container.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse.

IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Rinse mouth. Do NOT induce vomiting. Absorb spillage to prevent material damage.

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician.

STORAGE & DISPOSAL

Store locked up at room temperature. Store in original container or other corrosive resistant container. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com **Made in USA**

PK1273 Nail-RX RAPID AW-2026 Rev A 2020-04-29