



 **endure™** **Hematoxylin XT**  
H&E Staining System

REF #ED54

500ML CS/4

LOT



**For In Vitro Diagnostics Use**

86°F (30°C)  
59°F (15°C) 

Date Received: \_\_\_\_\_ Date Opened: \_\_\_\_\_

The Endure H&E System includes Hematoxylin XT ED5, Eosin ED2, Bluing ED3, and Clarifier ED4. Each component is analytically optimized for staining consistency and longevity.



**DANGER**  
**HAZARD STATEMENTS**

Causes skin irritation, and serious eye damage. May cause damage to organs (kidney) through prolonged or repeated exposure.

**PRECAUTIONARY STATEMENTS - PREVENTION**

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Do not breathe dust, vapors.

**PRECAUTIONARY STATEMENTS - RESPONSE**

**IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a doctor. Call a doctor if you feel unwell.

**IF ON SKIN:** Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

**STORAGE & DISPOSAL**

Store locked up at room temperature. Dispose of contents/container to an approved waste disposal plant.

**NOTICE:** Read Safety Data Sheet for additional information.

**EMERGENCY CONTACT:** Chemtrec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

[www.cancerdiagnostics.com](http://www.cancerdiagnostics.com)

**Made in USA**

P01492 Endure Hematoxylin XT 500mL Case AW-10412 Rev B 2020-09-29