



 **endure™** Hematoxylin XT
H&E Staining System

REF #ED5

500ML

LOT



For In Vitro Diagnostics Use

86°F (30°C)
59°F (15°C) 

Date Received: _____ Date Opened: _____

The Endure H&E System includes Hematoxylin XT ED5, Eosin ED2, Bluing ED3, and Clarifier ED4. Each component is analytically optimized for staining consistency and longevity.



DANGER
HAZARD STATEMENTS

Causes skin irritation, and serious eye damage. May cause damage to organs (kidney) through prolonged or repeated exposure.

PRECAUTIONARY STATEMENTS - PREVENTION

Wash body thoroughly after handling. Wear protective gloves. Wear eye and face protection. Do not breathe dust, vapors.

PRECAUTIONARY STATEMENTS - RESPONSE

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a doctor. Call a doctor if you feel unwell.

IF ON SKIN: Wash with plenty of water. Specific treatment, see first aid section on SDS. If skin irritation or rash occurs: get medical attention. Take off all contaminated clothing and wash it before reuse.

STORAGE & DISPOSAL

Store locked up at room temperature. Dispose of contents/container to an approved waste disposal plant.

NOTICE: Read Safety Data Sheet for additional information.

EMERGENCY CONTACT: Chemtec 1-800-424-9300



4300 Emperor Blvd - 400 Durham, NC 27703

Tel: 1-877-846-5393 Fax: 1-877-817-1716

www.cancerdiagnostics.com

Made in USA

PK1491 Endure Hematoxylin XT 500mL AW-10411 Rev.B 2020-09-29