

Nitrile Powder Free Examination Glove

Thanks to another advance in nitrile technology, Omni International is proud to announce its latest Nitrile Examination Glove, the OmniTrust™ Nitrile Glove. It's a new, thinner, lighter glove with impressive tactile sensitivity. You will find the OmniTrust™ Nitrile soft, stretchy and comfortable to wear with the comparable strength and protection of a regular nitrile glove. The glove maintains its natural durability for easy donning and improved elasticity. In addition, the textured fingers allow for excellent grip control and maximum ease of use with delicate instruments making the #201 glove the perfect choice for physicians, dentists and veterinarians alike. Experience the "OmniInnovation".



- Comfortable, Soft and Forms Like a Second Skin
- Not Made With Natural Rubber Latex
- Ambidextrous, Beaded Cuff, Single Use Only
- Textured Fingers for Tactile Sensitivity

Tensile Strength - Aged/Unaged	Elongation - Aged/Unaged
min. 14 MPa / min. 15 MPa	min. 400% / min. 500%

Dimensions	Finger	Palm	Cuff	Length	Width
Small	min. 0.07 mm 2.76 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	80 ± 10 mm
Medium	min. 0.07 mm 2.76 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	95 ± 10 mm
Large	min. 0.07 mm 2.76 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	110 ± 10 mm
Extra Large	min. 0.07 mm 2.76 mil	min. 0.06 mm 2.36 mil	min. 0.05 mm 1.97 mil	min. 240 mm	120 ± 10 mm

Item Numbers	Small	Medium	Large	X-Large
201 Series	201-121 201-221	201-122 201-222	201-123 201-223	201-124 201-224

Packaging	Boxes	Cases
packaged by weight	100 Gloves	10 Boxes
	200 Gloves*	10 Boxes

