


Schiff's Solution




Section 1. Identification

1.1 Product Identifier	Product Name: Schiff's Solution Product Number: SSC1117	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com

Section 2. Hazards Identification

2.1 GHS Classification	Skin irritation (Category 2) – H315 Serious eye damage (Category 1) – H318 Carcinogenicity (Category 1B) – H350	
Classification: Directive 1999/45/EC		C R22, R23, R31, R41, R45; S9, S26, S36/37/39, S45, S46 EN: Corrosive DA: Ætsende DE: Ätzend EL: Διαβρωτικό ES: Corrosivo FR: Corrosif IT: Corrosivo NL: Bijtend PT: Corrosivo SV: Frätande

Schiff's Solution

	 <p>Xn EN: Harmful DA: Sundhedsskadelig DE: Gesundheitsschädlich EL: Επιβλαβές ES: Nocivo FR: Nocif IT: Nocivo NL: Schadelijk PT: Nocivo SV: Hälsoskadlig</p> <p>For full text of R- & S- phrases: see Section 15.</p>								
2.2 Label Elements	<p>Pictogram(s):</p>  <p>Signal word: Danger</p> <p>Hazard statement(s): H318 – Causes serious eye damage H315 – Causes skin irritation H390 – May cause cancer H335 – May cause respiratory irritation.</p> <p>Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Get immediate medical advice/attention.</p>								
NFPA Scale: 0-4									
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr> <td>HEALTH</td> <td>2</td> </tr> <tr> <td>FLAMMABILITY</td> <td>0</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>G</td> </tr> </table>	HEALTH	2	FLAMMABILITY	0	PHYSICAL HAZARD	0	PERSONAL PROTECTION	G
HEALTH	2								
FLAMMABILITY	0								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	G								
2.3 Other Hazards	<p>PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.</p>								







Schiff's Solution

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Sodium Metabisulfite	7681-57-4	231-673-0	 Danger. 1 H318 Causes serious eye damage.  Warning. 4 H302 Harmful if swallowed.	≤ 2
Hydrochloric Acid	7647-01-0	231-595-7	 Warning. 3. H335 May cause respiratory irritation.  Danger. 1B. H314 Causes severe skin burns and eye damage.	≤ 1
Pararosaniline HCl (Basic Red 9)	569-61-9	209-321-2	 Warning. 2 H315 Causes skin irritation. 2 H319 Causes serious eye irritation. 3 H335 May cause respiratory irritation.  Danger. 1A H350 May cause cancer	≤ 0.5

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Schiff's Solution

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 2-8°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Sodium Metabisulfite: <i>NIOSH REL</i> : TWA 5 mg/m ³ OSHA PEL ‡: none Hydrochloric Acid: <i>NIOSH REL</i> : C 5 ppm (7 mg/m ³) OSHA PEL : C 5 ppm (7 mg/m ³)
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Schiff's Solution

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Transparent
Odor	Pungent Sulfur Dioxide
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Bases, oxidizers, sodium nitrate.
10.6 Hazardous Decomposition Materials	Oxides of sulfur.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Corrosive to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity.
---	---

Schiff's Solution

	No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). Hydrochloric Acid is listed as Group 3 – Not classifiable as to carcinogenicity in humans IARC: Pararosaniline (C.I. Basic red 9) is listed as a Group 2B: Possibly carcinogenic to humans. National Toxicology Program (NTP). Reasonably anticipated to be a human carcinogen.


Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA, IMDG, ADR	DOT: 3265 ADR/RID: 3265 IMDG: 3265 IATA: 3265
14.2 UN Proper Shipping Name DOT, IATA, IMDG, ADR	ADR/RID: CORROSIVE LIQUID, ACIDIC, ORGANIC, N.O.S. (Hydrochloric acid) IMDG: CORROSIVE LIQUID, ACIDIC, ORGANIC, N.O.S. (Hydrochloric acid) IATA: CORROSIVE LIQUID, ACIDIC, ORGANIC, N.O.S. (Hydrochloric acid) DOT: Corrosive liquid, acidic, organic, n.o.s. (Hydrochloric acid)
14.3 Transport Hazard Class(es)	DOT: 8 IATA: 8 IMDG: 8 ADR/RID: 8 
14.4 Packing Group	ADR/RID: III

Schiff's Solution

DOT, IATA,IMDG, ADR	IMDG: III IATA: III DOT: III
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.
14.7 Transport in Bulk According to Annex II of MARPOL73/78 and the IBC Code.	This product is not intended to be transported in bulk as defined in Annex II of MARPOL73/78 and the IBC Code.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Sara 313 Components	Hydrochloric Acid CAS# 7647-01-0
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	Pararosaniline: Warning! This product contains a chemical known to the State of California to cause cancer.
Right To Know Components	Pennsylvania, Massachusetts, New Jersey Right to Know: Sodium Metabisulphite CAS# 7681-57-4 Hydrochloric Acid CAS# 7647-01-0 Massachusetts, New Jersey Right to Know: Pararosaniline CAS# 569-61-9
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	<p>R22 EN: Harmful if swallowed. DA: Farlig ved indtagelse. DE: Gesundheitsschädlich beim Verschlucken. EL: Επιβλαβές σε περίπτωση καταπόσεως. ES: Nocivo por ingestión. FR: Nocif en cas d'ingestion. IT: Nocivo per ingestione. NL: Schadelijk bij opname door de mond. PT: Nocivo por ingestão. SV: Farligt vid förtäring.</p> <p>R23 EN: Toxic by inhalation. DA: Giftig ved indånding. DE: Giftig beim Einatmen. EL: Τοξικό όταν εισπνέεται. ES: Tóxico por inhalación. FR: Toxique par inhalation. IT: Tossico per inalazione. NL: Vergiftig bij inademing. PT: Tóxico por inalação. SV: Giftigt vid inandning.</p> <p>R31 EN: Contact with acids liberates toxic gas. DA: Udvikler giftig gas ved kontakt med syre. DE: Entwickelt bei Berührung mit Säure giftige Gase.</p>

Schiff's Solution

	<p>EL: Σε επαφή με οξέα ελευθερώνονται τοξικά αέρια. ES: Σε επαφή με οξέα ελευθερώνονται τοξικά αέρια. FR: Au contact d'un acide, dégage un gaz toxique. IT: A contatto con acidi libera gas tossico. NL: Vormt vergiftige gassen in contact met zuren. PT: Em contacto com ácidos liberta gases tóxicos. SV: Utvecklar giftig gas vid kontakt med syra.</p> <p>R41 EN: Risk of serious damage to eyes. DA: Risiko for alvorlig øjensskade. DE: Gefahr ernster Augenschäden. EL: Κίνδυνος σοβαρών οφθαλμικών βλαβών. ES: Riesgo de lesiones oculares graves. FR: Risque de lésions oculaires graves. IT: Rischio di gravi lesioni oculari. NL: Gevaar voor ernstig oogletsel. PT: Risco de lesões oculares graves. SV: Risk för allvarliga ögonskador.</p> <p>R45 EN: May cause cancer. DA: Kan fremkalde kræft. DE: Kann Krebs erzeugen. FR: Peut provoquer le cancer. IT: Può provocare il cancro. ES: Puede causar cáncer. FR: Risque de lésions oculaires graves. IT: Rischio di gravi lesioni oculari. NL: Gevaar voor ernstig oogletsel. PT: Risco de lesões oculares graves. SV: Risk för allvarliga ögonskador</p>
<p>S phrases</p>	<p>S9 EN: Keep container in a well-ventilated place. DA: Emballagen skal opbevares på et godt ventileret sted. DE: Behälter an einem gut gelüfteten Ort aufbewahren. EL: Το δοχείο να διατηρείται σε καλά αεριζόμενο μέρος. ES: Consérvese el recipiente en lugar bien ventilado. FR: Conserver le récipient dans un endroit bien ventilé. IT: Conservare il recipiente in luogo ben ventilato. NL: Op een goed geventileerde plaats bewaren. PT: Manter o recipiente num local bem ventilado. SV: Förpackningen förvaras på väl ventilerad plats.</p> <p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen. PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista. SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S36/37/39</p>

Schiff's Solution

	<p>EN: Wear suitable protective clothing, gloves and eye/face protection. DA: Brug særligt arbejdstøj, egnede beskyttelseshandsker og -briller/ansigtsskærm. DE: Bei der Arbeit geeignete Schutzkleidung, Schutzhandschuhe und Schutzbrille/Gesichtsschutz tragen. EL: Φοράτε κατάλληλη προστατευτική ενδυμασία, γάντια και συσκευή προστασίας ματιών/προσώπου ES: Usense indumentaria y guantes adecuados y protección para los ojos/la cara. FR: Porter un vêtement de protection approprié, des gants et un appareil de protection des yeux/du visage. IT: Usare indumenti protettivi e guanti adatti e proteggersi gli occhi/la faccia. NL: Draag geschikte beschermende kleding, handschoenen en een beschermingsmiddel voor de ogen/het gezicht. PT: Usar vestuário de protecção, luvas e equipamento protector para os olhos /face adequados. SV: Använd lämpliga skyddskläder, skyddshandskar samt skyddsglasögon eller ansiktsskydd.</p> <p>S45 EN: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). DA: Ved ulykkestilfælde eller ved ildebefindende er omgående lægebehandling nødvendig; vis etiketten, hvis det er muligt. DE: Bei Unfall oder Unwohlsein sofort Arzt zuziehen (wenn möglich, dieses Etikett vorzeigen). EL: Σε περίπτωση ατυχήματος ή αν αισθανθείτε αδιαθεσία ζητήστε αμέσως ιατρική συμβουλή (δείξτε την ετικέτα αν είναι δυνατό). ES: En caso de accidente o malestar, acúdase inmediatamente al médico (si es posible, muéstresele la etiqueta). FR: En cas d'accident ou de malaise, consulter immédiatement un medecin (si possible lui montrer l'etiquette). IT: In caso di incidente o di malessere consultare immediatamente il medico (se possibile, mostrargli l'etichetta). NL: Bij een ongeval of indien men zich onwel voelt, onmiddellijk een arts raadplegen (indien mogelijk hem dit etiket tonen). PT: Em caso de acidente ou de indisposição, consultar imediatamente o médico (se possível mostrar-lhe o rótulo). SV: Vid olycksfall, illamående eller annan påverkan, kontakta omedelbart läkare. Visa om möjligt etiketten.</p> <p>S46 EN: If swallowed, seek medical advice immediately and show this container or label. DA: Ved indtagelse, kontakt odømgående læge og vis denne beholder eller etiket. DE: Bei Verschlucken sofort ärztlichen Rat einholen und Verpackung odor Etikett vorzeigen. EL: Σε περίπτωση κατάποσης ζητήστε αμέσως ιατρική συμβουλή και δείξτε αυτό το δοχείο ή την ετικέτα. ES: En caso de ingestión, acúdase inmediatamente al medico y muéstresele la etiqueta o el envase. FR: En cas d'ingestion, consulter immédiatement un médecin et lui montrer l'emballage ou l'étiquette. IT: In caso d'ingestione consultare immediatamente il medico e mostrargli il contenitore o l'etichetta. NL: In geval van inslikken onmiddellijk een arts raadplegen en verpakking of etiket tonen. PT: Em caso de ingestão, consultar imediatamente o médico e mostrar-lhe a embalagem ou o rótulo. SV: Vid förtäring kontakta genast läkare och visa denna förpackning eller etiketten.</p>
--	--

Section 16. Other Information

Revision Date: 2019-01-11


Notice to Reader:



To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Periodic Acid Solution Product Number: SSC1103	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 Classification: Regulation (EC) No 1272/2008 (CLP)	Skin Irritation (Category 2), H315 Serious eye damage (Category 1), H318	
Classification: OSHA 29 CFR 1910.1200 HCS	Skin Irritation (Category 2), H315 Serious eye damage (Category 1), H318	
Classification: Directive 1999/45/EC		C, R38, R41, S26, S36/37/39 EN: Corrosive DA: Ætsende DE: Ätzend EL: Διαβρωτικό ES: Corrosivo FR: Corrosif IT: Corrosivo NL: Bijtend PT: Corrosivo SV: Frätande For full text of R- & S- phrases: see Section 15.


2.2 Label Elements	Pictogram(s):  Signal word: Warning Hazard statement(s): H315 – Causes skin irritation H318 – Causes serious eye damage Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water. P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes. P337 + 313 – If eye irritation persists get medical advice/attention.								
NFPA Scale: 0-4									
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr> <td>HEALTH</td> <td>1</td> </tr> <tr> <td>FLAMMABILITY</td> <td>0</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>B</td> </tr> </table>	HEALTH	1	FLAMMABILITY	0	PHYSICAL HAZARD	0	PERSONAL PROTECTION	B
HEALTH	1								
FLAMMABILITY	0								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	B								
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.								

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Water	7732-18-5	231-791-2	n/a	≥ 99
Periodic Acid	10450-60-9	233-937-0	 Danger. 2. H272 May intensify fire; oxidizer. Danger. 1B, H314 Causes severe skin burns and eye damage.	≤ 1

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention if irritation persists.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get

medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 2-8°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: N/A
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product.

	Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Colorless
Odor	Odorless
Odor Threshold	Unknown
pH	2.0 ± 0.5
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Oxidizing liquid

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Strong reducing agents, Strong bases, Powdered metals, Dimethyl sulfoxide (DMSO)
10.6 Hazardous Decomposition Materials	None known.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Corrosive to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA, IMDG, ADR	DOT: Not dangerous goods ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods
14.2 UN Proper Shipping Name DOT, IATA, IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods

14.3 Transport Hazard Class(es)	DOT: Not dangerous goods IATA: Not dangerous goods IMDG: Not dangerous goods ADR/RID: Not dangerous goods
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	<p>R38 EN: Irritating to skin. DA: Irriterer huden. DE: Reizt die Haut. EL: Ερεθίζει το δέρμα. ES: Irrita la piel. FR: Irritant pour la peau. IT: Irritante per la pelle. NL: Irriterend voor de huid. PT: Irriterend voor de huid. SV: Irriterar huden.</p> <p>R41 EN: Risk of serious damage to eyes. DA: Risiko for alvorlig øjenskaade. DE: Gefahr ernster Augenschäden. EL: Κίνδυνος σοβαρών οφθαλμικών βλαβών. ES: Riesgo de lesiones oculares graves. FR: Risque de lésions oculaires graves. IT: Rischio di gravi lesioni oculari. NL: Gevaar voor ernstig oogletsel. PT: Risco de lesões oculares graves. SV: Risk för allvarliga ögonskador.</p>
S phrases	<p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico.</p>

	<p>FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste.</p> <p>IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico.</p> <p>NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen.</p> <p>PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista.</p> <p>SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S36/37/39</p> <p>EN: Wear suitable protective clothing, gloves and eye/face protection.</p> <p>DA: Brug særligt arbejdstøj, egnede beskyttelseshandsker og -briller/ansigtsskærm.</p> <p>DE: Bei der Arbeit geeignete Schutzkleidung, Schutzhandschuhe und Schutzbrille/Gesichtsschutz tragen.</p> <p>EL: Φοράτε κατάλληλη προστατευτική ενδυμασία, γάντια και συσκευή προστασίας ματιών/προσώπου</p> <p>ES: Usense indumentaria y guantes adecuados y protección para los ojos/la cara.</p> <p>FR: Porter un vêtement de protection approprié, des gants et un appareil de protection des yeux/du visage.</p> <p>IT: Usare indumenti protettivi e guanti adatti e proteggersi gli occhi/la faccia.</p> <p>NL: Draag geschikte beschermende kleding, handschoenen en een beschermingsmiddel voor de ogen/het gezicht.</p> <p>PT: Usar vestuário de protecção, luvas e equipamento protector para os olhos /face adequados.</p> <p>SV: Använd lämpliga skyddskläder, skyddshandskar samt skyddsglasögon eller ansiktsskydd.</p>
--	---

Section 16. Other Information

Revision Date: 2019-01-11


Notice to Reader:



To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Bluing Reagent Product Number: SSC1152	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 GHS Classification	Skin irritation (Category 2) – H315 Eye irritation (Category 2) – H320	
Classification: Directive 1999/45/EC	 <p>Xi, R36/38; S26/39 EN: Irritant DA: Lokalirriterende DE: Reizend EL: Ερεθιστικό ES: Irritante FR: Irritant IT: Irritante NL: Irriterend PT: Irritante SV: Irriterande</p> <p>For full text of R- & S- phrases: see Section 15.</p>	




2.2 Label Elements	Pictogram(s):  Signal word: Warning Hazard statement(s): H315 – Causes skin irritation H320 – Causes eye irritation Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes.		
NFPA Scale: 0-4			
HMIS (U.S.A.) Scale: 0-4	HEALTH	1	
	FLAMMABILITY	0	
	PHYSICAL HAZARD	0	
	PERSONAL PROTECTION	B	
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.		

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Tris	77-86-1	201-064-4	 Warning. 2 H315 Causes skin irritation. 2 H319 Causes serious eye irritation. 3 H335 May cause respiratory irritation.	≤ 1
Sodium Hydroxide	1310-73-2	215-185-5	 Warning. 2 H315 Causes skin irritation. 2 H319 Causes serious eye irritation. 3 H335 May cause respiratory irritation.  Danger. 1A. H314 Causes severe skin burns and eye damage.	≤ 1

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Sodium Hydroxide: NIOSH REL: C 2 mg/m ³ OSHA PEL \pm : TWA 2 mg/m ³
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.
Engineering Controls	Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Clear
Odor	Odorless
Odor Threshold	Unknown
pH	9.0 \pm 0.10
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Metals, oxidizers, reducing agents, acids, alkalis.
10.6 Hazardous Decomposition Materials	None known.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Irritating to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods. ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods.
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods. IATA: Not dangerous goods. IMDG: Not dangerous goods. ADR/RID: Not dangerous goods.
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	None of the components in this mixture are listed.
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R36/38 EN: Irritating to eyes and skin. DA: Irriterer øjnene og huden DE: Reizt die Augen und die Haut. EL: Ερεθίζει τα μάτια και το δέρμα. ES: Irrita los ojos y la piel. FR: Irritant pour les yeux et la peau. IT: Irritante per gli occhi e la pelle. NL: Irriterend voor de ogen en de huid. PT: Irritante para os olhos e pele. SV: Irriterar ögonen och huden.
S phrases	S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e

	<p>consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen. PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista. SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S37/39 EN: Wear suitable gloves and eye/face protection. DA: Brug egnede beskyttelseshandsker og -briller/ansigtsskærm under arbejdet. DE: Bei der Arbeit geeignete Schutzhandschuhe und Schutzbrille/Gesichtsschutz tragen. EL: Φοράτε κατάλληλα γάντια και συσκευή προστασίας ματιών/προσώπου. ES: Úsense guantes adecuados y protección para los ojos/la cara. FR: Porter des gants appropriés et un appareil de protection des yeux/du visage. IT: Usare guanti adatti e proteggersi gli occhi/la faccia. NL: Draag geschikte handschoenen en een beschermingsmiddel voor de ogen/het gezicht. PT: Usar luvas e equipamento protector para os olhos /face adequados. SV: Använd lämpliga skyddshandskar samt skyddsglasögon eller ansiktsskydd.</p>
--	---

Section 16. Other Information

Revision Date: 2019-01-11


Notice to Reader:



To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Alpha-Amylase Solution (1%) Product Number: SSC1058 Product Synonym: Diastase Solution	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	877-846-5393
	Fax Number	877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com

Section 2. Hazards Identification

2.1 GHS Classification	Respiratory sensitisation (Category 1) – H334 Acute toxicity, Oral (Category 4) – H302	
Classification: Directive 1999/45/EC	 Xn R20/22, R42; S24, S28, S36/37, S45, S61 EN: Harmful DA: Sundhedsskadelig DE: Gesundheitsschädlich EL: Επιβλαβές ES: Nocivo FR: Nocif IT: Nocivo NL: Schadelijk PT: Nocivo SV: Hälsoskadlig For full text of R- & S- phrases: see Section 15.	




2.2 Label Elements	Pictogram(s):  Signal word: Danger Hazard statement(s): H334 – May cause allergy or asthma symptoms or breathing difficulties if inhaled. H303 – May be harmful if swallowed. Precautionary statement(s): P261 – Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray. P342 + P311 – If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.								
NFPA Scale: 0-4									
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr> <td>HEALTH</td> <td>1</td> </tr> <tr> <td>FLAMMABILITY</td> <td>0</td> </tr> <tr> <td>PHYSICAL HAZARD</td> <td>0</td> </tr> <tr> <td>PERSONAL PROTECTION</td> <td>B</td> </tr> </table>	HEALTH	1	FLAMMABILITY	0	PHYSICAL HAZARD	0	PERSONAL PROTECTION	B
HEALTH	1								
FLAMMABILITY	0								
PHYSICAL HAZARD	0								
PERSONAL PROTECTION	B								
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.								

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Sodium Azide	26628-22-8	247-852-1	 Danger 1, 2 H300 Fatal if swallowed., H310 Fatal in contact with skin.  Warning. 1. H410 Very toxic to aquatic life with long lasting effects.	≤ 0.1
Sodium Chloride	7647-14-5	231-598-3	 Warning. 1. H317 May cause an allergic skin reaction.	≤ 0.8
Alpha-Amylase	9000-90-2	232-565-6	n/a	≤ 1.0

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed
 See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed
 No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosive hazards expected
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures
 Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions
 Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up
 Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.
 Avoid contact with skin and eyes.
 Wash thoroughly after handling.
 Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.
 Store in well ventilated area.
 Keep container tightly closed.
 Store at 2-8°C.

7.3 Specific end use(s).
 See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: Sodium Azide ACGIH TLV: 0.29/mg/m ³ ceiling NIOSH REL-C: 0.3mg/m ³ skin
8.2 Exposure controls	Personal Protective Equipment (PPE): Eye/Face protection.

	<p>Safety glasses or goggles are required.</p> <p>Skin protection. Protective clothing is required.</p> <p>Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove.</p> <p>Respiratory protection. Avoid breathing vapor.</p> <p>Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
Engineering Controls	<p>Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.</p>

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Colorless to slightly yellow.
Odor	Odorless
Odor Threshold	Unknown
pH	Unknown
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	Sodium Azide may react with lead and copper plumbing to form explosive metal azides.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Strong acids, bases and oxidizers.
10.6 Hazardous Decomposition Materials	None known.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Irritating to eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity Sodium Azide is a known mutagen.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA, IMDG, ADR	DOT: Not dangerous goods ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods
14.2 UN Proper Shipping Name DOT, IATA, IMDG, ADR	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods IATA: Not dangerous goods IMDG: Not dangerous goods ADR/RID: Not dangerous goods
14.4 Packing Group	ADR/RID: Not dangerous goods

DOT, IATA,IMDG, ADR	IMDG: Not dangerous goods IATA: Not dangerous goods DOT: Not dangerous goods
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	Sodium Azide is listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	Sodium Azide Massachusetts Pennsylvania New Jersey Sodium Chloride Pennsylvania New Jersey Alpha-Amylase Pennsylvania New Jersey
Other Regulations	Sodium Azide SARA 302 SARA 313
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R20/22 EN: Harmful by inhalation and if swallowed. DA: Farlig ved indånding og ved indtagelse. DE: Gesundheitsschädlich beim Einatmen und Verschlucken. EL: Επιβλαβές όταν εισπνέεται και σε περίπτωση καταπόσεως. ES: Nocivo por inhalación y por ingestión. FR: Nocif par inhalation et par ingestion. IT: Nocivo per inalazione e ingestione. NL: Schadelijk bij inademing en opname door de mond. PT: Nocivo por inalação e ingestão. SV: Farligt vid inandning och förtäring. R42 EN: May cause sensitization by inhalation. DA: Kan give overfølsomhed ved indånding. DE: Sensibilisierung durch Einatmen möglich. EL: Μπορεί να προκαλέσει ευαισθητοποίηση όταν εισπνέεται. ES: Posibilidad de sensibilización por inhalación. FR: Peut entraîner une sensibilisation par inhalation. IT: Può provocare sensibilizzazione per inalazione. NL: Kan overgevoeligheid veroorzaken bij inademing. PT: Pode causar sensibilização por inalação. SV: Kan ge allergi vid inandning.
S phrases	S24 EN: Avoid contact with skin. DA: Undgå kontakt med huden. DE: Berührung mit der Haut vermeiden.

EL: Αποφεύγετε την επαφή με το δέρμα.
 ES: Evítese el contacto con la piel.
 FR: Éviter le contact avec la peau.
 IT: Evitare il contatto con la pelle.
 NL: Aanraking met de huid vermijden.
 PT: Evitar o contacto com a pele.
 SV: Undvik kontakt med huden.

S28

EN: After contact with skin, wash immediately with plenty of water.
 DA: Kommer stof på huden vaskes straks med store mængder ... (angives af fabrikanten).
 DE: Bei Berührung mit der Haut sofort abwaschen mit viel ... (vom Hersteller anzugeben).
 EL: Σε περίπτωση επαφής με το δέρμα, πλυθείτε αμέσως με άφθονο ... (το είδος του υγρού καθορίζεται από τον παραγωγό).
 ES: En caso de contacto con la piel, lávese inmediata y abundantemente con ... (productos a especificar por el fabricante).
 FR: Après contact avec la peau, se laver immédiatement et abondamment avec ... (produits appropriés à indiquer par le fabricant).
 IT: In caso di contatto con la pelle lavarsi immediatamente ed abbondantemente con ... (prodotti idonei da indicarsi da parte del fabbricante).
 NL: Na aanraking met de huid onmiddellijk wassen met veel ... (aan te geven door de fabrikant).
 PT: Após contacto com a pele, lavar imediata e abundantemente com ... (produtos adequados a indicar pelo produtor).
 SV: Vid kontakt med huden tvätta genast med mycket ... (anges av tillverkaren).

S36/37

EN: Wear suitable protective clothing and gloves.
 DA: Brug særligt arbejdstøj og egnede beskytteshandsker.
 DE: Bei der Arbeit geeignete Schutzhandschuhe und Schutzkleidung tragen.
 EL: Φοράτε κατάλληλη προστατευτική ενδυμασία και γάντια.
 ES: Usense indumentaria y guantes de protección adecuados.
 FR: Porter un vêtement de protection et des gants appropriés.
 IT: Usare indumenti protettivi e guanti adatti.
 NL: Draag geschikte handschoenen en beschermende kleding.
 PT: Usar vestuário de protecção e luvas adequadas.
 SV: Använd lämpliga skyddskläder och skyddshandskar.

S45

EN: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
 DA: Ved ulykkestilfælde eller ved ildebefindende er omgående lægebehandling nødvendig; vis etiketten, hvis det er muligt.
 DE: Bei Unfall oder Unwohlsein sofort Arzt zuziehen (wenn möglich, dieses Etikett vorzeigen).
 EL: Σε περίπτωση ατυχήματος ή αν αισθανθείτε αδιαθεσία ζητήστε αμέσως ιατρική συμβουλή (δείξτε την ετικέτα αν είναι δυνατό).
 ES: En caso de accidente o malestar, acúdase inmediatamente al médico (si es posible, muéstresele la etiqueta).
 FR: En cas d'accident ou de malaise, consulter immédiatement un medecin (si possible lui montrer l'etiquette).
 IT: In caso di incidente o di malessere consultare immediatamente il medico (se possibile, mostrargli l'etichetta).
 NL: Bij een ongeval of indien men zich onwel voelt, onmiddellijk een arts raadplegen (indien mogelijk hem dit etiket tonen).
 PT: Em caso de acidente ou de indisposição, consultar imediatamente o médico (se possível mostrar-lhe o rótulo).
 SV: Vid olycksfall, illamående eller annan påverkan, kontakta omedelbart läkare. Visa om möjligt etiketten.

S61

EN: Avoid release to the environment. Refer to special instructions/Safety data sheets.
 DA: Undgå udledning til miljøet. Se særlig vejledning/leverandørbrugsanvisning.
 DE: Freisetzung in die Umwelt vermeiden. Besondere Anweisungen einholen/Sicherheitsdatenblatt zu Rate ziehen.
 EL: Αποφύγετε την ελευθέρωσή του στο περιβάλλον. Αναφερθείτε σε ειδικές οδηγίες/.ελτίο δεδομένων ασφαλείας.
 ES: Evítese su liberación al medio ambiente. Recábense instrucciones específicas de la ficha de datos de seguridad.
 FR: Éviter le rejet dans l'environnement. Consulter les instructions spéciales/la fiche de données de sécurité.

	<p>IT: Non disperdere nell'ambiente. Riferirsi alle istruzioni speciali/ schede informative in materia di sicurezza.</p> <p>NL: Voorkom lozing in het milieu. Vraag om speciale instructies/veiligheidskaart.</p> <p>PT: Evitar a libertação para o ambiente. Obter instruções específicas/fichas de segurança.</p> <p>SV: Undvik utsläpp till miljön. Läs särskilda instruktioner/varuinformationsblad.</p>
--	--

Section 16. Other Information

Revision Date: 2019-01-11


Notice to Reader:



To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.

Section 1. Identification

1.1 Product Identifier	Product Name: Hematoxylin, Mayer's (Lillie's Modification) Product Number: SSC1059	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com/sds

Section 2. Hazards Identification

2.1 GHS classification	Acute toxicity, Oral (Category 4) – H302 Skin irritation (Category 2) – H315 Eye irritation (Category 2A) – H319
Classification: Directive 1999/45/EC	 <p>Xi, R36/38; S2425, S26, S45 EN: Irritant DA: Lokalirriterende DE: Reizend EL: Ερεθιστικό ES: Irritante FR: Irritant IT: Irritante NL: Irriterend PT: Irritante SV: Irriterande</p> <p>For full text of R- & S- phrases: see Section 15.</p>





2.2 Label Elements	Pictogram(s):  Signal word: Warning Hazard statement(s): H302 – Harmful if swallowed. H315 – Causes skin irritation H320 – Causes serious eye irritation Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 - IF IN EYES: Rinse cautiously with water for several minutes. P301 + P312 - IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.										
NFPA Scale: 0-4											
HMIS (U.S.A.) Scale: 0-4	<table border="1"> <tr><td>HEALTH</td><td>1</td></tr> <tr><td>FLAMMABILITY</td><td>0</td></tr> <tr><td>PHYSICAL HAZARD</td><td>0</td></tr> <tr><td>PERSONAL PROTECTION</td><td>B</td></tr> </table>	HEALTH	1	FLAMMABILITY	0	PHYSICAL HAZARD	0	PERSONAL PROTECTION	B		
HEALTH	1										
FLAMMABILITY	0										
PHYSICAL HAZARD	0										
PERSONAL PROTECTION	B										
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.										



Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Glycerin	56-81-5	200-289-5	n/a	≤ 30
Aluminum Sulfate	7784-31-8	n/a	 Warning. 2 H315 Causes skin irritation. 3 H335 May cause respiratory irritation.  Danger. 2 H318 Causes serious eye damage.  Warning. 2 H411 Toxic to aquatic life with long lasting effects.	≤ 5
Acetic Acid	64-19-7	200-580-7	 Warning. 3, H226	≤ 2

			 Flammable liquid and vapour. Danger. 1A H314 Causes severe skin burns and eye damage.	
Hematoxylin	517-28-2	208-237-3	 Warning. 4 H302 Harmful if swallowed. 2 H315 Causes skin irritation. 2 H319 Causes serious eye irritation. 3 H335 May cause respiratory irritation.	≤ 1

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage



<p>7.1 Precautions for safe handling. Avoid contact with skin and eyes. Wash thoroughly after handling. Avoid breathing vapor.</p> <p>7.2 Conditions for safe storage, including any incompatibilities. Store in well ventilated area. Keep container tightly closed. Store at 15-30°C.</p> <p>7.3 Specific end use(s). See section 1.2</p>

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	<p>Exposure Limits: Acetic acid, glacial; ACGIH: 15ppm STEL (37 mg/m³ STEL) 10 ppm TWA (25 mg/m³ TWA); NIOSH: 15 ppm STEL (37 mg/m³ STEL) 10 ppm TWA (25 mg/m³ TWA); OSHA: 10 ppm TWA (25 mg/m³ TWA)</p>
8.2 Exposure controls	<p>Personal Protective Equipment (PPE): Eye/Face protection. Safety glasses or goggles are required. Skin protection. Protective clothing is required. Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove. Respiratory protection. Avoid breathing vapor. Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
Engineering Controls	<p>Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.</p>

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Purple
Odor	Odorless
Odor Threshold	Unknown
pH	2.4 ± 0.3
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown

Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-Ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	No relevant data available.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.
10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Reactive with strong oxidizing agents, reducing agents, metals, acids, and alkalis.
10.6 Hazardous Decomposition Materials	Carbon Monoxide. Carbon Dioxide. May release irritating or toxic fumes under fire conditions.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data available Skin Corrosion/Irritation. Irritating to skin and mucous membranes. Serious Eye Damage/Irritation. Irritating to eyes. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data available.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: Acetic Acid: LD50 Oral, Rat, 3310 mg/kg. LD50 Dermal, Rabbit, 1060 mg/kg.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.

12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA,IMDG, ADR	DOT: Not dangerous goods. ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods.
14.2 UN Proper Shipping Name DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.3 Transport Hazard Class(es)	DOT: Not dangerous goods. IATA: Not dangerous goods. IMDG: Not dangerous goods. ADR/RID: Not dangerous goods.
14.4 Packing Group DOT, IATA,IMDG, ADR	ADR/RID: Not dangerous goods. IMDG: Not dangerous goods. IATA: Not dangerous goods. DOT: Not dangerous goods.
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.
14.7 Transport in Bulk According to Annex II of MARPOL73/78 and the IBC Code.	This product is not intended to be transported in bulk as defined in Annex II of MARPOL73/78 and the IBC Code.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous	None of the components in this mixture are listed.

Substances; Section 355	
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	This mixture does not contain a chemical known to the state of California to cause cancer.
Right To Know Components	<p>Acetic Acid: California Massachusetts Minnesota New Jersey Pennsylvania Rhode Island</p> <p>Glycerin: Illinois Massachusetts Minnesota Pennsylvania Rhode Island</p> <p>Acetic Acid is listed on Canada's Ingredient Disclosure List.</p>
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	<p>R36/38 EN: Irritating to eyes and skin. DA: Irriterer øjnene og huden DE: Reizt die Augen und die Haut. EL: Ερεθίζει τα μάτια και το δέρμα. ES: Irrita los ojos y la piel. FR: Irritant pour les yeux et la peau. IT: Irritante per gli occhi e la pelle. NL: Irriterend voor de ogen en de huid. PT: Irritante para os olhos e pele. SV: Irriterar ögonen och huden.</p>
S phrases	<p>S24/25 EN: Avoid contact with skin and eyes. DA: Undgå kontakt med huden og øjnene. DE: Berührung mit den Augen und der Haut vermeiden. EL: Αποφεύγετε επαφή με το δέρμα και με τα μάτια. ES: Evítese el contacto con los ojos y la piel. FR: Éviter le contact avec la peau et les yeux. IT: Evitare il contatto con gli occhi e con la pelle. NL: Aanraking met de ogen en de huid vermijden. PT: Evitar o contacto com a pele e os olhos. SV: Undvik kontakt med huden och ögonen.</p> <p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávense inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen.</p>

	<p>PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista.</p> <p>SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S45</p> <p>EN: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).</p> <p>DA: Ved ulykkestilfælde eller ved ildebefindende er omgående lægebehandling nødvendig; vis etiketten, hvis det er muligt.</p> <p>DE: Bei Unfall oder Unwohlsein sofort Arzt zuziehen (wenn möglich, dieses Etikett vorzeigen).</p> <p>EL: Σε περίπτωση ατυχήματος ή αν αισθανθείτε αδιαθεσία ζητήστε αμέσως ιατρική συμβουλή (δείξτε την ετικέτα αν είναι δυνατό).</p> <p>ES: En caso de accidente o malestar, acúdase inmediatamente al médico (si es posible, muéstresele la etiqueta).</p> <p>FR: En cas d'accident ou de malaise, consulter immédiatement un medecin (si possible lui montrer l'etiquette).</p> <p>IT: In caso di incidente o di malessere consultare immediatamente il medico (se possibile, mostrargli l'etichetta).</p> <p>NL: Bij een ongeval of indien men zich onwel voelt, onmiddellijk een arts raadplegen (indien mogelijk hem dit etiket tonen).</p> <p>PT: Em caso de acidente ou de indisposição, consultar imediatamente o médico (se possível mostrar-lhe o rótulo).</p> <p>SV: Vid olycksfall, illamående eller annan påverkan, kontakta omedelbart läkare. Visa om möjligt etiketten.</p>
--	--

Section 16. Other Information

Revision Date: 2019-01-11

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.