

250 ml Histology Container

Product Code: YC250



General Description


- Ideal for specimens and general purpose storage and transportation
- Shatter-resistant, translucent body
- Screw on lid creates secure seal

Material:

Polypropylene Body, FDA Certified
HDPE Cap, FDA Certified
Product is recyclable

Packing Details

Case:	100 pc/case
Packaged:	10/sleeve, 10 sleeves
Case Dimensions:	19 x 8 x 16"
Case Weight:	10.75 lbs.
Cube Dimension:	1.41
Cases per Skid:	48



Dimensions listed are in millimeters