

# 8 oz Pathology Container (Tall)

**Product Code:** PC8250-T

## General Description

Color: Translucent  
Body Diameter: 84.2 mm  
Body Height: 68.9 mm




- Snap-on lid creates secure seal
- Shatter-resistant, translucent body allows for easy viewing

### **Material:**

HDPE, FDA Certified  
Country of Origin: China  
*Product is recyclable*



Case: 250 sets/case  
Packaged: 250 Bagged  
Case Dimensions: 16.50 x 6.50 x 15.25"  
Case Weight: 9.8 lbs  
Cube Dimension: 0.95  
Cases per Skid: 60



Product dimensions listed are in millimeters