



CRYOSTAT SPECIMEN CHUCKS

**Fast-Freeze core, with deep-grooved  
OCT “grab patterns”**

O-ring colors help eliminate sample errors.  
Orientation mark assures proper sample positioning.

### Cryostar



**99-3530F** Thermo  
Flat Circular 30mm



**99-3530FW** Thermo  
Flat Waffle 30mm



**99-3540F** Thermo  
Flat Circular 40mm



**99-3540FW** Thermo  
Flat Waffle 40mm

### Avantik QS/Microm HM Series



**SC0030** Avantik Stem  
Circular-30mm



**SC0030W** Avantik Stem  
Waffle-30mm



**SC0040** Avantik Stem  
Circular-40mm



**SC0040W** Avantik Stem  
Waffle-40mm

### Leica (except for CMI1950) TBS, Tanner



**99-3322** Leica Stem  
Circular-22mm



**99-3325** Leica Stem  
Circular-25mm



**99-3330** Leica Stem  
Circular-30mm



**99-3330W** Leica  
Stem Waffle-30mm



**99-3340** Leica Stem  
Circular-40mm



**99-3340W** Leica Stem  
Waffle-40mm



**99-3355** Leica Stem  
Circular-55mm



### Leica CMI1950



**99-3430F** Leica Flat  
Circular 30mm



**99-3430FW** Leica Flat  
Waffle 30mm



**99-3440F** Leica Flat  
Circular 40mm



**99-3440FW** Leica Flat  
Waffle 40mm

**Shandon**



**99-3630** Shandon Stem  
Circular-30mm




**99-3630W** Shandon Stem  
Waffle-30mm







**99-3640** Shandon Stem  
Circular-40mm



**99-3640W** Shandon Stem  
Waffle-40mm



**FROZEN SECTION  
SPECIMEN CHUCKS**

**CryoCore™ Circular and Waffle Tissue Chucks**  
 Cryostat embedding chucks are made of anodized aluminum with a circular or waffle pattern for secure tissue mounting. Color coded for multiple specimens. Specify desired O-Ring color by adding (-1,-2,-3,-4,-5,-6,-7,-8,-9 after item #). White, Blue, Teal, Red, Yellow, Orange, Purple, Light Pink and Black (-1, -2, -3, -4, -5, -6, -7, -8, -9)

Description	Size	Waffle	Circular	Colors	White	Blue	Teal	Red	Yellow	Orange	Purple	Light Pink	Black
					1	2	3	4	5	6	7	8	9
Avantik Stem	30mm	SC0030W	SC0030	1-8	○	●	●	●	●	●	●	●	●
Avantik Stem	40mm	SC0040W	SC0040	1-9	○	●	●	●	●	●	●	●	●
Leica Stem	22mm		99-3322	1-8	○	●	●	●	●	●	●	●	●
Leica Stem	25mm		99-3325	1-8	○	●	●	●	●	●	●	●	●
Leica Stem	30mm	99-3330W	99-3330	1-8	○	●	●	●	●	●	●	●	●
Leica Stem	40mm	99-3340W	99-3340	1-9	○	●	●	●	●	●	●	●	●
Leica Stem	55mm		99-3355	1-8	○	●	●	●	●	●	●	●	●
Leica Flat	30mm	99-3430FW	99-3430F	1-8	○	●	●	●	●	●	●	●	●
Leica Flat	40mm	99-3440FW	99-3440F	1-9	○	●	●	●	●	●	●	●	●
Shandon Stem	30mm	99-3630W	99-3630	1-8	○	●	●	●	●	●	●	●	●
Shandon Stem	40mm	99-3640W	99-3640	1-9	○	●	●	●	●	●	●	●	●
Thermo Flat	30mm	99-3530FW	99-3530F	1-8	○	●	●	●	●	●	●	●	●
Thermo Flat	40mm	99-3540FW	99-3540F	1-9	○	●	●	●	●	●	●	●	●