



Nitric Acid Solution 1%



Section 1. Identification

1.1 Product Identifier	Product Name: Nitric Acid Solution (1%) Product Number: SSC1201	
1.2 Intended use	EN: Laboratory reagent. For professional use only. DA: Laboratoriereagens. Kun til professionelt brug. DE: Laboratoriumreagens. Alleen voor professioneel gebruik. EL: Αντιδραστήριο εργαστηρίου. Για επαγγελματική χρήση μόνο. ES: Reactivo de laboratorio. Sólo para uso profesional. FR: Réactif de laboratoire. Pour un usage professionnel uniquement. IT: Laboratorio di reagente. Solo per uso professionale. NL: Laboratoriumreagens. Alleen voor professioneel gebruik. PT: Reagente de laboratório. Para uso profissional. SV: Laboratoriereagens. Endast för yrkesmässig användning.	
1.3 Details of the supplier of the safety data sheet	Manufacturer	Cancer Diagnostics, Inc.
	Address	116 Page Point Circle Durham, NC 27703
	Phone Number	1-877-846-5393
	Fax Number	1-877-817-1716
	e-mail	compliance@cancerdiagnostics.com
	Website	www.cancerdiagnostics.com

Section 2. Hazards Identification

2.1 GHS Classification	Oxidizing liquids (Category 3) – H272 Skin corrosion (Category 1A) – H314 Serious eye damage (Category 1) – H 318	
Classification: Directive 1999/45/EC	 C R34, R41; S25, S26 S36/37/39 EN: Corrosive DA: Ætsende DE: Ätzend EL: Διαβρωτικό ES: Corrosivo FR: Corrosif IT: Corrosivo NL: Bijtend PT: Corrosivo SV: Frätande	 O:

Nitric Acid Solution 1%


	EN: Oxidizing DA: Brandnærende DE: Brandfördernd EL: Οξειδωτικό ES: Comburente FR: Comburant IT: Comburente NL: Oxyderend PT: Comburente SV: Oxiderande For full text of R- & S- phrases: see Section 15.		
2.2 Label Elements	Pictogram(s):  Signal word: Danger Hazard statement(s): H314 – Causes severe skin burns and eye damage H272 – May intensify fire; oxidizer Precautionary statement(s): P302+ P352 - IF ON SKIN: Wash with plenty of soap and water P305 + P351 + P335 - IF IN EYES: Rinse cautiously with water for several minutes. Get immediate medical advice/attention.		
NFPA Scale: 0-4			
HMIS (U.S.A.) Scale: 0-4	HEALTH	2	
	FLAMMABILITY	0	
	PHYSICAL HAZARD	1	
	PERSONAL PROTECTION	B	
2.3 Other Hazards	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.		

Section 3. Composition and Information on Ingredients

3.2 Chemical Description: Mixture

*May contain additional non-hazardous proprietary ingredients.

*May contain additional active ingredients at concentrations <0.1%w/v.

Hazardous Ingredients:	CAS#	EC#	GHS Symbols	%
Nitric Acid	7697-37-2	231-714-2	 Danger. 3 H272 May intensify fire; oxidizer. 1 H314 Causes severe skin burns and eye damage. 1 H318 Causes serious eye damage.	≤ 1.0

Nitric Acid Solution 1%

Section 4. First Aid Measures

4.1 Description of first aid measures

Eye Contact: Check for and remove contact lenses. Immediately flush eyes with copious amounts of water and get immediate medical attention.

Skin Contact: Remove contaminated clothing and wash contact area with mild soap and copious amounts of water. Get medical attention if irritation develops.

Inhalation: If inhaled, remove person to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention if symptoms worsen.

Ingestion: Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Loosen tight clothing such as collar, tie, belt or waistband. Get immediate medical attention.

4.2 Most important symptoms and effects, both acute and delayed

See section 2.2 and 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available

Section 5. Fire Fighting Measures

5.1 Extinguishing Media	Extinguish fire using water spray, carbon dioxide, chemical foam, or dry chemical.
5.2 Special hazards arising from the substance or mixture	No unusual fire or explosion hazards expected.
5.3 Advice for firefighters	As with any fire, wear personal protection equipment, including a self-contained breathing apparatus (S.C.B.A.)

Section 6. Accidental Release Measures

6.1 Personal precautions, protective equipment and emergency procedures

Wear chemical resistant clothing, gloves, and eye protection. Wear NIOSH/MSHA approved breathing apparatus.

6.2 Environmental precautions

Keep material away from heat, flame, ignition sources, and reactive materials. Don't allow product to enter drain.

6.3 Methods and materials for containment and cleaning up

Wipe up or absorb spill using inert absorbent and place in a suitable waste container for disposal.

Section 7. Handling and Storage

7.1 Precautions for safe handling.

Avoid contact with skin and eyes.

Wash thoroughly after handling.

Avoid breathing vapor.

7.2 Conditions for safe storage, including any incompatibilities.

Store in well ventilated area.

Keep container tightly closed.

Store at 15-30°C.

7.3 Specific end use(s).

See section 1.2

Section 8. Exposure Controls / Personal Protection

8.1 Control parameters	Exposure Limits: N/A
-------------------------------	--------------------------------

Nitric Acid Solution 1%

8.2 Exposure controls	<p>Personal Protective Equipment (PPE):</p> <p>Eye/Face protection. Safety glasses or goggles are required.</p> <p>Skin protection. Protective clothing is required.</p> <p>Hand protection. Chemical resistant gloves are required. Glove material must be resistant to the components of this product. Consult glove manufacturer for specific recommendations of appropriate material and thickness of glove.</p> <p>Respiratory protection. Avoid breathing vapor.</p> <p>Environmental exposure controls. Avoid releasing large quantities into the environment. No additional information.</p>
Engineering Controls	<p>Working area should be adequately large and well ventilated to prevent concentration of vapors. Provide mechanical exhaust ventilation or other engineering controls to keep airborne concentrations of vapors below their respective threshold limits.</p>

Section 9. Physical and Chemical Properties

Physical State	Liquid
Color	Colorless
Odor	Odorless
Odor Threshold	Unknown
pH	Acidic
Melting Point/ Freezing Point	Unknown
Initial Boiling Point	Unknown
Flash Point	Unknown
Evaporation Rate	Unknown
Flammability (solid, gas)	Unknown
Upper/Lower Flammability Limits	Unknown
Vapor Pressure	Unknown
Vapor Density	Unknown
Relative Density	Unknown
Solubility(ies)	Water
Partition Coefficient: n-octanol/water	Unknown
Auto-ignition Temperature	Unknown
Decomposition Temperature	Unknown
Viscosity	Unknown
Explosive Properties	Not explosive.
Oxidizing Properties	Unknown

Section 10. Stability and Reactivity

10.1 Reactivity	Stable under normal temperatures and pressures.
10.2 Chemical Stability	Stable under normal temperatures and pressures.
10.3 Possibility of Hazardous Reactions	No hazardous reactions known.

Nitric Acid Solution 1%

10.4 Conditions to Avoid	Fire, static electricity, direct sunlight.
10.5 Incompatible Materials	Strong acids, bases and oxidizers
10.6 Hazardous Decomposition Materials	None known.

Section 11. Toxicological Information

11.1 Information on Toxicological Effects.	Acute Toxicity. No relevant data Skin Corrosion/Irritation. Corrosive to skin and mucous membranes. Serious Eye Damage/Irritation. Corrosive to the eye. Respiratory or skin sensitization. No relevant data available. Germ Cell Mutagenicity. No relevant data.
Carcinogenicity.	International Agency for Research on Cancer (IARC). None of the components are listed. National Toxicology Program (NTP). None of the components are listed.

Section 12. Ecological Information

12.1 Toxicity	Fish: No relevant studies identified. Crustacea: No relevant studies identified. Algae/Aquatic Plants: No relevant studies identified. Other Organisms: No relevant studies identified.
12.2 Persistence and Degradability.	No relevant studies identified.
12.3 Bioaccumulative Potential.	No relevant studies identified.
12.4 Mobility in Soil.	Miscible in water. May spread in water systems. This component is non-volatile.
Additional Remarks	None.
12.5 Results of PBT and vPvB Assessment.	PBT: This mixture does not contain any substances that are assessed to be a PBT. vPvB: This mixture does not contain any substances that are assessed to be a vPvB.

Section 13. Disposal Considerations

13.1 Waste Disposal Methods.	Sewage disposal is discouraged. Waste should not be disposed of by release to sewers. Dispose waste in accordance with federal, state and local environmental control regulations.
Product/Packaging Disposal.	Final decisions on the appropriate waste management method must be in line with local, regional and national regulations.
Other Disposal Recommendations.	No relevant data available.

Section 14. Transport Information

14.1 UN Number DOT, IATA, IMDG, ADR	DOT: 2031 ADR/RID: 2031 IMDG: 2031 IATA: 2031
--	--

Nitric Acid Solution 1%

14.2 UN Proper Shipping Name DOT, IATA, IMDG, ADR	ADR/RID: NITRIC ACID IMDG: NITRIC ACID IATA: NITRIC ACID DOT: NITRIC ACID
14.3 Transport Hazard Class(es)	DOT: 8 IATA: 8 IMDG: 8 ADR/RID: 8
14.4 Packing Group DOT, IATA, IMDG, ADR	ADR/RID: II IMDG: II IATA: II DOT: II
14.5 Environmental Hazards	Marine Pollutant: No
14.6 Special Precautions for User	Not applicable.

Section 15. Regulatory Information

15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture.	
Extremely Hazardous Substances; Section 355	None of the components in this mixture are listed.
Toxic Substances Control Act; TSCA	All of the components in this mixture are listed.
California Proposition 65	None of the components in this mixture are listed.
Right To Know Components	<u>Nitric Acid</u> Massachusetts Pennsylvania New Jersey
15.2 Chemical Safety Assessment	No Chemical Safety Assessment has been carried out for this substance/mixture by Cancer Diagnostics, Inc.
R phrases	R34 EN: Causes burns. DA: Ætsningsfare. DE: Verursacht Verätzungen. EL: Προκαλεί εγκαύματα. ES: Provoca quemaduras. FR: Provoque des brulures. IT: Provoca ustioni. NL: Veroorzaakt brandwonden. PT: Provoca queimaduras. SV: Frätande. R41 EN: Risk of serious damage to eyes. DA: Risiko for alvorlig øjenskade. DE: Gefahr ernster Augenschäden. EL: Κίνδυνος σοβαρών οφθαλμικών βλαβών. ES: Riesgo de lesiones oculares graves. FR: Risque de lésions oculaires graves. IT: Rischio di gravi lesioni oculari. NL: Gevaar voor ernstig oogletsel. PT: Risco de lesões oculares graves. SV: Risk för allvarliga ögonskador.
S phrases	S25 EN: Avoid contact with eyes. DA: Undgå kontakt med øjnene.

Nitric Acid Solution 1%

	<p>DE: Berührung mit den Augen vermeiden. EL: Αποφεύγετε την επαφή με τα μάτια. ES: Evítese el contacto con los ojos. FR: Éviter le contact avec les yeux. IT: Evitare il contatto con gli occhi. NL: Aanraking met de ogen vermijden. PT: Evitar o contacto com os olhos. SV: Undvik kontakt med ögonen.</p> <p>S26 EN: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. DA: Kommer stoffet i øjnene, skylles straks grundigt med vand og læge kontaktes. DE: Bei Berührung mit den Augen sofort gründlich mit Wasser abspülen und Arzt konsultieren. EL: Σε περίπτωση επαφής με τα μάτια πλύνετε τα αμέσως με άφθονο νερό και ζητήστε ιατρική συμβουλή. ES: En caso de contacto con los ojos, lávensse inmediata y abundantemente con agua y acúdase a un médico. FR: En cas de contact avec les yeux, laver immédiatement et abondamment avec de l'eau et consulter un spécialiste. IT: In caso di contatto con gli occhi, lavare immediatamente e abbondantemente con acqua e consultare un medico. NL: Bij aanraking met de ogen onmiddellijk met overvloedig water afspoelen en deskundig medisch advies inwinnen. PT: Em caso de contacto com os olhos, lavar imediata e abundantemente com água e consultar um especialista. SV: id kontakt med ögonen, spola genast med mycket vatten och kontakta läkare.</p> <p>S36/37/39 EN: Wear suitable protective clothing, gloves and eye/face protection. DA: Brug særligt arbejdstøj, egnede beskyttelseshandsker og -briller/ansigtsskærm. DE: Bei der Arbeit geeignete Schutzkleidung, Schutzhandschuhe und Schutzbrille/Gesichtsschutz tragen. EL: Φοράτε κατάλληλη προστατευτική ενδυμασία, γάντια και συσκευή προστασίας ματιών/προσώπου ES: Úsense indumentaria y guantes adecuados y protección para los ojos/la cara. FR: Porter un vêtement de protection approprié, des gants et un appareil de protection des yeux/du visage. IT: Usare indumenti protettivi e guanti adatti e proteggersi gli occhi/la faccia. NL: Draag geschikte beschermende kleding, handschoenen en een beschermingsmiddel voor de ogen/het gezicht. PT: Usar vestuário de protecção, luvas e equipamento protector para os olhos /face adequados. SV: Använd lämpliga skyddskläder, skyddshandskar samt skyddsglasögon eller ansiktsskydd.</p>
--	---

Section 16. Other Information

Revision Date: 2019-01-11

Notice to Reader:

To the best of our knowledge, the information contained herein is accurate. However, neither the above named supplier nor any of its subsidiaries assumes any liability whatsoever for the accuracy or completeness of the information contained herein.