



4x
MAGNIFICATION
CAMERA

PATHCO
a Cancer Diagnostics, Inc. division

EDG
ENHANCED DIAGNOSTIC GROSSING



See what's new
Point your camera at this
QR code

EDG GS5 Grossing Station

equipment@cancerdiagnostics.com


cancerdiagnostics.com



EDG

ENHANCED DIAGNOSTIC GROSSING

Precision Airflow. Seamless Construction. A Better Grossing Station.



- PATHCO Air-Curtain Ventilation
- Seamless Formed Stainless Steel Construction
- Optimized Airflow Containment

Engineered airflow designed to help limit formalin vapor migration. A unique airflow barrier positioned at the operator interface helps contain vapors and directs airflow downward across the work surface toward the rear exhaust system.

Choose from different built-in imaging system options

PATHCO Imaging

Easy to configure
Easy to deliver
Easy to install



Founded in 1960 by a pathologist, PATHCO introduced instruments and blades designed for the grossing bench and autopsy suite. Today, PATHCO continues that tradition with grossing stations engineered for the modern pathology laboratory.



PATHCO, a Division of Cancer Diagnostics Inc.