



# ESPO SLIDE PRINTER

ON DEMAND PRINTER FOR TISSUE SLIDES

**With the launch of the ESPO Path Slide Printer, an on demand slide printer for anatomical pathology labs, Cancer Diagnostics, Inc. sets a new unmatched milestone in slide printing technology.**

Despite its minimal size, the ESPO Path Slide Printer has two independent magazines for slides to automatically print on two types of glass. A front loading feeder is also available to print on separate slides already with tissue. This allows positive immune control to be labeled directly on the slide.

#### FEATURES:

**“Risk management” and “workability” are compatible**

- Front loading for single slides
- Ultra-compact size
- Touch panel display
- WiFi, LAN, USB connection
- Excellent print quality

#### TECHNICAL DATA:

<b>Print method</b>	Line Thermal transfer method
<b>Print dot density</b>	300dpi
<b>Printing speed</b>	20 slides/min
<b>Print direction</b>	0°, 90°, 180°, 270°
<b>Print quality</b>	Water - Alcohol - Xylene tolerant
<b>Barcode type</b>	UPC-E, Code 39, Code 93, Code 128, EAN 128, NW-7, QR Code, Micro QR Code, DMC
<b>Supply of slide</b>	3-way supply, 2 x magazine, 1 x front loading
<b>Slide removal</b>	Stacker for exclusive order of each one (20 sheets)
<b>Frost width</b>	15mm, 20mm, 25mm
<b>Ink Ribbon</b>	Exclusive solvent resistant ink ribbon
<b>Interface</b>	Wired LAN, Wireless LAN (IEEE802.11b/g/n), USB 2.0



**Matsunami Japan** (Established 1844) ISO 13485

  
 Dimensions  
 254mm (Height)  
 241mm (Long)  
 127mm (wide)

---

  
 Weight  
 ~6.35kg



<b>Cancer Diagnostics Inc.</b>	<b>Contact Information</b>	<b>Visit us online</b>
--------------------------------	----------------------------	------------------------

4300 Emperor Blvd-400  
Durham, North Carolina 27703

**Toll Free:** (877) 846 5393  
**Fax:** (877) 817 1716

[cancerdiagnostics.com](http://cancerdiagnostics.com)