

SafeCube Dispenser™

**Reduce The Risk Of Accidental Drips, Drops & Spills!
& Enhances Regulatory Compliance!**

The Safest method for filling smaller containers with laboratory chemicals from a five gallon cubitainer, eg., 10% Formalin, Disinfectants, & Reagents.

System Benefits

- No more Dripping Spigots all over the floor from a 5-gallon cube hanging over the edge of the counter.
- No need to balance a container while filling, simply place container to be filled in filling tray and turn on spigot.
- When not filling containers close it up for safe storage
- No more spills from hooking spigots with lab coats.
- Control the Drips, Drops & Accidental Spills with Secondary Spill Containment.
(Main Cabinet holds up to five (5) Gallons)
- Liquid Detection Alarm immediately signals working personnel of any accidental spills.

Perfect For Use In:

- Surgery and Labor & Delivery
- Ambulatory Surgery
- Histo-Pathology Laboratory
- Clinical and Research Laboratory
- College & University Laboratory
- Morgue & Autopsy Suites
- Housekeeping & Central Services

PRODUCT # SCD-100A
Measurements: 20" deep x 18 1/2" wide x 17" high



American Bio-Safety, LLC

Bus; (916) 652-8021 4322 Anthony Court, # 5 Rocklin, CA. 95677 Fax; (916) 652-8020

WebSite; www.AmericanBioSafety.com Toll Free; 800-624-8021 Email; absafety@sbcglobal.net

Cubitainer® Combination Package and Cube® Insert are Registered Trademarks of the Hedwin Corporation, Baltimore, Maryland.