

# SafeCube Dispenser™

**Reduce The Risk Of Accidental Drips, Drops & Spills!  
& Enhances Regulatory Compliance!**

**The Safest method for filling smaller containers with laboratory chemicals from a five gallon cubitainer, eg., 10% Formalin, Disinfectants, & Reagents.**

## System Benefits

- No more Dripping Spigots all over the floor from a 5-gallon cube hanging over the edge of the counter.
- No need to balance a container while filling, simply place container to be filled in filling tray and turn on spigot.
- When not filling containers close it up for safe storage
- No more spills from hooking spigots with lab coats.
- Control the Drips, Drops & Accidental Spills with Secondary Spill Containment.  
(Main Cabinet holds up to five (5) Gallons)
- Liquid Detection Alarm immediately signals working personnel of any accidental spills.

## Perfect For Use In:

- Surgery and Labor & Delivery
- Ambulatory Surgery
- Histo-Pathology Laboratory
- Clinical and Research Laboratory
- College & University Laboratory
- Morgue & Autopsy Suites
- Housekeeping & Central Services

PRODUCT # SCD-1000  
Measurements: 20" deep x 18 ½" wide x 17" high

