

# 500 ml Histology Container

**Product Code:** YC500



## General Description

- Ideal for specimens and general purpose storage and transportation
- Shatter-resistant, translucent body
- Screw on lid creates secure seal


## **Material:**

Polypropylene Body, FDA Certified  
HDPE Cap, FDA Certified

*Product is recyclable*

## Packing Details

Case:	100 pc/case
Packaged:	10/sleeve, 10 sleeves
Case Dimensions:	19 x 8 x 18.25"
Case Weight:	12.9 lbs.
Cube Dimension:	1.61
Cases per Skid:	48



Dimensions listed are in millimeters