






















SYMBOL	SYMBOL TITLE	EXPLANATORY TEXT
	Manufacturer	Symbol indicates the medical device manufacturer.
	Date of manufacture	Symbol indicates the date when the medical device was manufactured.
	Authorized European representative	Symbol indicates the Authorized representative in the European Community.
	Catalogue or model number	Symbol indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial number	Symbol indicates the manufacturer's serial number so that a specific medical device can be identified.
	Batch code	Symbol indicates the manufacturer's batch code so that the batch or lot can be identified.
	Use by	Symbol indicates the date after which the medical device is not to be used.
	In vitro diagnostic medical device	Symbol indicates a medical device that is intended to be used as an in vitro diagnostic medical device.
	Consult instructions for use	Symbol indicates the need for the user to consult the instructions for use.
	Caution: Read all warnings and precautions in instructions for use	Symbol indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Storage temperature range	Symbol indicates the temperature limits to which the medical device can be safely exposed.
	Keep dry	Symbol indicates a medical device that needs to be protected from moisture.

SYMBOL	SYMBOL TITLE	EXPLANATORY TEXT
	Fragile, handle with care	Symbol indicates a medical device that can be broken or damaged if not handled carefully.
	Do not reuse	Symbol indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Do not resterilize	Symbol indicates a medical device that is not to be resterilized.
	Contains natural rubber or latex	Symbol indicates the presence of natural rubber or dry natural rubber latex as a material of construction within the medical device or the packaging of a medical device.
	Non sterile	Symbol indicates a medical device that has not been subjected to a sterilization process. No symbol. Always represented in text.
	Product subjected to a sterilization process	Symbol indicates a medical device that has been subjected to a sterilization process.
	Sterilized by ethylene oxide treatment	Symbol indicates a medical device that has been sterilized using ethylene oxide.
	Sterilized using irradiation	Symbol indicates a medical device that has been sterilized using irradiation.
	Sterilized using steam or dry heat	Symbol indicates a medical device that has been sterilized using steam or dry heat.
	Recycle: Electronic Equipment	DO NOT THROW IN TRASH.